

SIGVARIS
GROUP

Compression Solutions



Graduated Compression Socks,
Hosiery and Inelastic Garments

Established in 1864, SIGVARIS GROUP is, and has always been, a family-owned company. Founded and headquartered in Switzerland, our passion for helping people feel their best is why we've been a global market leader for 158 years and counting.

Our compression garment portfolio includes Medical, Consumer Health Care and products for Medically Complex Edema that are specifically designed to restore, assist and enhance people's lives. We are renowned for making exceptional products that combine advanced technology and quality craftsmanship with diverse styles that fit our modern customers' lives.

When heritage, competence, and dedication come together, there's no limit to where we'll go.

Live your life. Be yourself.

For more information, please visit:
sigvaris.com

TO ORDER:
800.363.4999 (phone)
800.263.8736 (fax)
canadaorders@sigvaris.com

Online Dealer Portal:
sigvaris-online.com

The World-Wide Leader in Compression Therapy

WHY SIGVARIS GROUP?

We produce the best compression garments to enhance our customers' health and quality of life. We are dedicated to helping legs look and feel their best. Using premium yarns, we've spun ourselves into being the world-wide leader in graduated compression socks and hosiery by embracing our core beliefs of QUALITY, COMFORT, FIT and INNOVATION. Staying true to these values, we are continuing our efforts by expanding our business to help customers with medically-complex needs.



COMFORT

We offer the most comfortable fabric options available, including the finest cottons, wools and microfibers. We double-cover the yarns for ease of donning and doffing.

QUALITY

More than 20 people handle our products from start to finish, ensuring every garment is of the highest quality.



PRECISE FIT

Our products are designed using the most recent human data measurement studies, fitting 99.49% of the population.

NATURAL FIBERS

We have the highest use of natural materials, including Sea Island Cotton, Supima Cotton, Linen and Merino Wool.



VARIETY

We offer you more choices of materials, colours, patterns, and styles than any other medical compression brand.

PANTYHOSE

Constructed with body-shaping contouring and a cotton lined gusset for breathability and comfort.



GRIP TOPS

High stretch, high memory components to deliver comfortable "stay up" power.

RECIPROCATED TOE

Many of our products are constructed with a smooth, flat toe seam and a contoured toe box for added comfort.

INNOVATION

We use the finest textiles and materials to create modern styles in our circular-knit products, while maintaining medical efficacy. We continue to invest in medical research to help people feel their best, every day.

CUSTOMER SERVICE

Canada-wide salesforce and a live customer care team combine to provide dedicated support & expertise.

	PRODUCT	STYLES	COMPRESSION LEVEL	PAGE
ESSENTIAL	COTTON 230	● calf, calf w/grip-top, thigh w/grip-top ○ calf, thigh w/grip-top	20–30 mmHg 30–40 mmHg	10
	OPAQUE 860	● calf, petite calf, calf plus, calf w/grip-top, thigh w/grip-top, petite thigh w/grip-top, extra-large thigh w/grip-top, pantyhose, petite pantyhose, plus-sized pantyhose ○ calf, calf plus, calf w/grip-top, thigh w/grip-top, extra-large thigh w/grip-top, thigh w/waist attachment	20–30 mmHg 30–40 mmHg	12, 14
STYLE	SHEER 780 (for women only)	● calf, thigh w/grip-top, pantyhose, maternity pantyhose ○ calf, thigh w/grip-top, pantyhose	15–20 mmHg 20–30 mmHg 30–40 mmHg	18
	MEDIUM SHEER 750 (for women only)	● calf, thigh w/grip-top, pantyhose	20–30 mmHg	20
	SOFT OPAQUE 840 (for women only)	● calf, calf w/grip-top, thigh w/grip-top, pantyhose, maternity pantyhose ○ calf, thigh w/grip-top, pantyhose	15–20 mmHg 20–30 mmHg 30–40 mmHg	22
	PATTERNS 710 (for women only)	● thigh w/grip-top, pantyhose	20–30 mmHg	24
	SEA ISLAND COTTON 220	● calf	20–30 mmHg	26
	MERINO WOOL 240	● calf	20–30 mmHg	28
	LINEN 250	● calf	20–30 mmHg	30
	MICROFIBER 820 (for men only)	● calf, calf w/grip-top, thigh w/grip-top, extended foot sizes	15–20 mmHg 20–30 mmHg 30–40 mmHg	32
	MICROFIBER PATTERNS 830	● calf	20–30 mmHg	34
	CUSHIONED COTTON 360	● calf	20–30 mmHg	38
MOTION	HIGH TECH 412C	● calf	20–30 mmHg	40
	THERMOREGULATING WOOL 422	● calf	20–30 mmHg	42
	COMFORT 850	● calf	20–30 mmHg	44
	NATURAL RUBBER 500	○ calf, thigh, thigh w/grip-top, thigh w/waist attachment, pantyhose	30–40 mmHg 40–50 mmHg 50–60 mmHg	48
SPECIALTY	ADVANCE ARMSLEEVE 910	○ armsleeve w/ grip-top w/gauntlet and without gauntlet	15–20 mmHg 20–30 mmHg 30–40 mmHg	50

Closed Toe = ● Open Toe = ○

	PRODUCT	STYLES	COMPRESSION LEVEL	PAGE
CONSUMER HEALTH	CASUAL COTTON 146/186	● calf	15–20 mmHg	54
	TRAVENO TRAVEL SOCKS 146/T - 186/T	● calf	15–20 mmHg	56
	CUSHIONED COTTON 142/182	● calf	15–20 mmHg	58
	SEA ISLAND COTTON 151/191	● calf	15–20 mmHg	60
	ALL-SEASON MERINO WOOL 152/192	● calf	15–20 mmHg	62
	LINEN 153/193	● calf	15–20 mmHg	64
	MICROFIBER SHADES 143/183	● calf	15–20 mmHg	66
	SHEER FASHION 120	● calf, thigh w/grip-top, pantyhose ○ calf	15–20 mmHg	68
	SOFT SILHOUETTE LEGGINGS 170	○ open foot	15–20 mmHg	70
	ATHLETIC RECOVERY SOCKS 401	● calf	15–20 mmHg	72
	PERFORMANCE SLEEVES 412V	○ open foot	15–20 mmHg	74
	MERINO OUTDOOR SOCKS 421	● calf	15–20 mmHg	76
	EVERSOFT DIABETIC SOCKS 160	● calf	8–15 mmHg	78
	MEDICALLY COMPLEX EDEMA & INELASTIC WRAPS			80
	ACCESSORIES			98
OTHER	SYMPTOMS & INDICATIONS			103
	GLOSSARY			104
	HOW TO MEASURE & ORDER			106

Natural fibers

Sigvaris uses premium natural fibers, including Sea Island Cotton, Merino Wool, and our latest addition, Linen. Ideal for sensitive skin, the natural fibers are soft, breathable, have a high moisture absorption capacity and are easy to care for.

sigvaris



Be yourself.

ESSENTIAL

Support Yourself

Sigvaris Essential solutions always meet your customers' needs. These straightforward stockings are an everyday go-to for people who seek the comfort, support, and reliability of traditional therapeutic socks and stockings.

- Why choose **Essential**?
- Reduce swelling and discomfort
 - Long-distance travel
 - Everyday comfort and reliability

STYLE

Express Yourself

Choice is empowering, and fashion can be as unique as you are. With Sigvaris Style, personality will shine as your customers select from distinctive hosiery and socks.

- Why choose **Style**?
- Fashionable, modern patterns and designs
 - Special occasion and formal wear
 - Our most luxurious feel

MOTION

Push Yourself

When daily comfort for hard-working legs is required, put your customers in motion. Sigvaris Motion socks offer all day strong support for people with active lives.

- Why choose **Motion**?
- Cushioned sole for added comfort
 - Durable thermoregulating wool
 - Moisture-wicking yarns





Support Yourself

Why choose **Essential**?

- Reduce swelling and discomfort
- Long-distance travel
- Everyday comfort and reliability

ESSENTIAL PRODUCTS

Series	Product	Page
230	Cotton	10
860	Opaque (Women)	12
860	Opaque (Men)	14

BEST SUITED FOR

- Containment of various types of edema, including venous edema, phlebolymphe~~de~~ma, and early stages of lymphedema
- Ideal for dress or casual shoes
- Suited for post-vein treatment



COTTON 230

for women & men

Compression Level

20–30 mmHg

30–40 mmHg

Features and Benefits



Sensitive Skin



Cotton



Breathable

- 100% Supima Cotton against the skin
- Ideal for sensitive skin and allergies
- Breathable and durable
- Excellent containment of various types of edema
- Classic ribbed design for everyday wear

Fiber Content

- 50% Nylon, 25% Spandex, 25% Supima Cotton (content can vary slightly by style)
- Latex free



20–30 mmHg
Closed Toe Calf w/ Grip-top
Medium Long in Black

232CMLM99/S



Silicone dotted grip-top available
for calf and thigh



Scan here to see product images online
sigvaris.com

SIZING CHART

	Small		Medium		Large		X Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")		29–36 cm (11.5"–14")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")		41–61 cm (16"–24")	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Thigh Circumference	45.5–66 cm (18"–26")		48–73 cm (19"–28.5")		52–84 cm (20.5"–33")		NA	
Leg Length: Women & Open Toe	Short	Long	Short	Long	Short	Long		
	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")		
Leg Length: Men	Short	Long	Short	Long	Short	Long		
	≤76 cm (30")	>76 cm (30")	≤76 cm (30")	>76 cm (30")	≤76 cm (30")	>76 cm (30")		
Size	SS	SL	MS	ML	LS	LL	XS	XL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12
Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS

Closed Toe = ● Open Toe = ○					
	WHITE 00	NAVY 10	BLACK MIST 14	LIGHT BEIGE 66	BLACK 99
Calf (C)*	●	●	●	●●	●
Calf w/Grip-Top (C/S) Silicone dotted grip-top band				●**	●
Thigh w/Grip-Top (N) Silicone dotted grip-top band				●●	●

* 30–40 mmHg extra large calf only available in black with grip-top
** 20–30 mmHg only

Product Code Example: 232NMLW66

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour
23	2	N	ML	W	66

Product Code Example: 232CLSM66/S

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour	Grip-top
23	2	C	LS	M	66	/S



OPAQUE 860

for women

Compression Level

20–30 mmHg	30–40 mmHg
------------	------------

Features and Benefits



Durable



Largest Size Range



Opaque

- Largest size range, from petite to plus
- Durable construction combined with a smooth opaque look
- Comfortable, stretchable material easily glides on and off
- Containment of various types of edema
- Plus size pantyhose are an alternative to maternity pantyhose

Fiber Content

- 64% Nylon, 36% Spandex (content can vary slightly by style)
- Latex free



Scan here to see product images online
sigvaris.com



20–30 mmHg
Closed Toe Calf Medium Long
in Black

862CMLW99



Silicone dotted grip-top available
for calf and thigh

SIZING CHART

	Small			Medium			Large		X Large	
Ankle Circumference	18–21.5 cm (7"–8.5")			21.5–25.5 cm (8.5"–10")			25.5–29 cm (10"–11.5")		29–36 cm (11.5"–14")	
Calf Circumference	28–39 cm (11"–15.5")			31.5–45 cm (12.5"–17.5")			35.5–51 cm (14"–20")		41–61 cm (16"–24")	
Calf Circumference: Calf Plus w/Grip-Top	NA			NA	39.5–53 cm (15.5"–21")		45.5–61 cm (18"–24")		51–71 cm (20"–28")	
Calf Length	Petite	Short	Long	Petite	Short	Long	Short	Long	Short	Long
	≤35.5 cm (14")	≤41 cm (16")	>41 cm (16")	≤35.5 cm (14")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Thigh Circumference	45.5–66 cm (18"–26")			48–73 cm (19"–28.5")			52–84 cm (20.5"–33")		64–95 cm (25"–37.5")	
Leg Length	Petite	Short	Long	Petite	Short	Long	Short	Long	Short	Long
	≤65 cm (25.5")	≤72 cm (28.5")	>72 cm (28.5")	≤65 cm (25.5")	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")
Hip Circumference	83–112.5 cm (32.5"–44.5")			87–127 cm (34.5"–50")			98–150.5 cm (38.5"–59.5")		NA	
Plus or Maternity Size Hip Circumference	≤142 cm (56")			≤157.5 cm (62")			≤178 cm (70")		NA	
Size	SP	SS	SL	MP	MS	ML	LS	LL	XS	XL
Calf Plus w/Grip-Top					MS/P	ML/P	LS/P	LL/P	XS/P	XL/P

Women's Shoe Size : Petite Leg 4–8, Short Leg 5.5–9.5, Long Leg 8–12

STYLES AND COLOURS

Closed Toe = ●
Open Toe = ○

	HONEY 33	GOLDEN 36	LIGHT BEIGE 66	BLACK 99
Calf (C)	●	●	○	●●
Calf Petite Sizes (C) Only available in 20–30 mmHg, not available with dotted grip-top band	●	●		●
Calf w/Grip-Top (C/S) Silicone dotted grip-top band			●○	●
Calf Plus w/Grip-Top (C/P) Only available with silicone dotted grip-top band			●○	
Thigh w/Grip-Top (N)** Silicone dotted grip-top band	●	●	○	●
Thigh Petite Sizes (N) Only available in 20–30 mmHg	●	●		●
Pantyhose (P)	●	●		●
Pantyhose Petite (P) Only available in 20–30 mmHg	●	●		●
Pantyhose Plus and Maternity (P/P)			●	●
Thigh w/Waist Attachment (W)* Indicate leg LEFT (/L) or RIGHT (/R)			○	

* Not available in XL with waist attachment

** Thigh w/Grip-Top (N) and Open Toe Thigh w/Grip-Top (N) sizes XS and XL are only available in Light Beige 66

Product Code Example: 863PMSW33

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour
86	3	P	MS	W	33

Product Code Example: 862CLLO66/P

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour	Grip-top with Plus size
86	2	C	LL	O	66	/P



OPAQUE 860

for men

Compression Level

- 20–30 mmHg
- 30–40 mmHg

Features and Benefits



Durable



Largest Size Range



Opaque

- Durable construction with a ribbed design on closed toe calf products
- Smooth opaque design on open-toe styles.
- Comfortable, stretchable material easily glides on and off
- Containment of various types of edema

Fiber Content

- 64% Nylon, 36% Spandex
(content can vary slightly by style)
- Latex free



20–30 mmHg
Open Toe Thigh Medium Long
in Light Beige

862NMLO66



Silicone dotted grip-top available
for calf and thigh



Scan here to see product images online
sigvaris.com

SIZING CHART

	Small		Medium		Large		X Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")		29–36 cm (11.5"–14")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")		41–61 cm (16"–24")	
Calf Circumference: Calf Plus w/Grip-Top	NA		39.5–53 cm (15.5"–21")		45.5–61 cm (18"–24")		51–71 cm (20"–28")	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Thigh Circumference	45.5–66 cm (18"–26")		48–73 cm (19"–28.5")		52–84 cm (20.5"–33")		64–95 cm (25"–37.5")	
Leg Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")
Size	SS	SL	MS	ML	LS	LL	XS	XL
Calf Plus w/Grip-Top			MS/P	ML/P	LS/P	LL/P	XS/P	XL/P

Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS

	RIBBED LIGHT BEIGE 66	RIBBED BLACK 99	SMOOTH LIGHT BEIGE 66	SMOOTH BLACK 99
Closed Toe = ● Open Toe = ○				
Calf (C)	●*	●*	○	○
Calf w/Grip-Top (C/S) Silicone dotted grip-top band			○	
Calf Plus w/Grip-Top (C/P) Only available with silicone dotted grip-top band			○	
Open Toe Thigh w/Grip-Top (N) Silicone dotted grip-top band			○	
Thigh w/Waist Attachment (W)** Indicate leg LEFT (/L) or RIGHT (/R)			○	

* 20–30mmHg only
** Not available in XL with waist attachment

Product Code Example: 862CMLM99

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour
86	2	C	ML	M	99

Product Code Example: 862CLLO66/P

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour	Grip-top with Plus size
86	2	C	LL	O	66	/P



Express Yourself

Why choose **Style**?

- Fashionable, modern patterns and designs
- Special occasion and formal wear
- Our most luxurious feel

STYLE PRODUCTS

Series	Product	Page
780	Sheer	18
750	Medium Sheer	20
840	Soft Opaque	22
710	Patterns	24
220	Sea Island Cotton	26
240	Merino Wool	28
250	Linen	30
820	Microfiber	32
830	Microfiber Patterns	34

BEST SUITED FOR

- Management of Chronic Venous Disease (CVD)
- Pre- and post-vein treatment with minimal edema
- Pregnancy
- People who prefer fashionable options
- Ideal for business or dress



SHEER 780

for women

Compression Level

- 15–20 mmHg
- 20–30 mmHg
- 30–40 mmHg

Features and Benefits



Sheer



Elegant



Breathable

- Fashionable for professional and dress wear
- Consistent sheerness throughout the leg
- Excellent breathability and moisture management
- Maternity pantyhose options available

Fiber Content

- 15–20 and 20–30 mmHg: 69% Nylon, 31% Spandex (content can vary slightly by style)
- 30–40 mmHg: 60% Nylon, 40% Spandex (content can vary slightly by style)
- Latex free



15–20 mmHg
Closed Toe Thigh Medium Long in
Dark Navy

781NMLW08



Sensinnov grip-top on thigh



Scan here to see product images online
sigvaris.com

SIZING CHART

	Small		Medium		Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")	
Calf Length	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Thigh Circumference	45.5–66 cm (18"–26")		48–73 cm (19"–28.5")		52–84 cm (20.5"–33")	
Leg Length	Short	Long	Short	Long	Short	Long
	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")
Hip Circumference	83–112.5 cm (32.5"–44.5")		87–127 cm (34.5"–50")		98–150.5 cm (38.5"–59.5")	
Maternity hip circumference	≤142 cm (56")		≤157.5 cm (62")		≤178 cm (70")	
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

STYLES AND COLOURS

	DARK NAVY 08	HONEY 33	WARM SAND 34*	GOLDEN 36	TOASTED ALMOND 37*	CAFÉ 73*	MOCHA 85*	NIGHT-SHADE 94	BLACK 99
Calf (C)	●	●●	●●	●●	●●	●●	●●	●	●●
Thigh w/Grip-Top (N) Sensinnov grip-top band	●	●●	●●	●●	●●	●●	●●	●	●●
Pantyhose (P)	●	●●	●●	●●	●●	●●	●●	●	●●
Maternity Pantyhose (M)*		●							●

* 20–30 mmHg only

Product Code Example: 781NMLW36

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour
78	1	N	ML	W	36

Product Code Example: 783PSSO33

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour
78	3	P	SS	O	33



MEDIUM SHEER 750

for women

Compression Level

20–30 mmHg

Features and Benefits



Medium Sheer



Durable



Breathable

- Fashionable, durable, medium sheer hosiery for everyday wear
- Excellent breathability and moisture management
- Comfortable 2" top band on calf style

Fiber Content

- 68% Nylon, 32% Spandex
(content can vary slightly by style)
- Latex free



20–30 mmHg
Closed Toe Thigh Medium Long in Black

752NMLW99



Silicone dotted grip-top on thigh



Scan here to see product images online
sigvaris.com

SIZING CHART

	Small		Medium		Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")	
Calf Length	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Thigh Circumference	45.5–66 cm (18"–26")		48–73 cm (19"–28.5")		52–84 cm (20.5"–33")	
Leg Length	Short	Long	Short	Long	Short	Long
	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")
Hip Circumference	83–112.5 cm (32.5"–44.5")		87–127 cm (34.5"–50")		98–150.5 cm (38.5"–59.5")	
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

STYLES AND COLOURS

Closed Toe = ●	<div><div>HONEY 33</div><div>NIGHT-SHADE 94</div><div>BLACK 99</div></div>		
	●	●	●
Calf (C)	●	●	●
Thigh w/Grip-Top (N) Silicone dotted grip-top band	●	●	●
Pantyhose (P)	●	●	●

Product Code Example: 752NMLW33

Series	mmHg Level	Style	Size	Gender	Colour
75	2	N	ML	W	33



SOFT OPAQUE 840

for women

Compression Level

- 15–20 mmHg
- 20–30 mmHg
- 30–40 mmHg

Features and Benefits



Opaque



Microfiber



Soft

- Fashionable opaque hosiery for professional and everyday wear
- Long-lasting microfiber softness provides all day comfort and allows product to easily glide on and off
- Durable and breathable
- Maternity pantyhose option available


Fiber Content

- 79% Nylon, 21% Spandex (content can vary slightly by style)
- Latex free



20–30 mmHg
Closed Toe Pantyhose Small Long in Black

842PSLW99



SensinnoV grip-top on thigh and silicone dotted grip-top available on calf



Scan here to see product images online
sigvaris.com

SIZING CHART

	Small		Medium		Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")	
Calf Length	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Thigh Circumference	45.5–66 cm (18"–26")		48–73 cm (19"–28.5")		52–84 cm (20.5"–33")	
Leg Length	Short	Long	Short	Long	Short	Long
	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")
Hip Circumference	83–112.5 cm (32.5"–44.5")		87–127 cm (34.5"–50")		98–150.5 cm (38.5"–59.5")	
Maternity hip circumference	≤142 cm (56")		≤157.5 cm (62")		≤178 cm (70")	
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

STYLES AND COLOURS

	MIDNIGHT BLUE 09	CHAI 35	PECAN 40	MULBERRY 60*	ESPRESSO 89*	GRAPHITE 91	BLACK 99
Closed Toe = ● Open Toe = ○							
Calf (C)	●	●●	●	●	●	●	●○
Calf w/Grip-Top (C/S)* Silicone dotted grip-top							●
Thigh w/Grip-Top (N) SensinnoV® grip-top band	●	●●	●	●	●	●	●○
Pantyhose (P)	●	●●	●	●	●	●	●○
Maternity Pantyhose (M)*							●

* 20–30 mmHg only

Product Code Example: 842CMSW09

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour
84	2	C	MS	W	09

PATTERNS 710

for women

Compression Level

20–30 mmHg

Features and Benefits



Sheer



Elegant



Breathable

- Stylish diamond-dot and stripe patterns for professional and dress wear
- Sheer hosiery with excellent breathability and moisture management
- Unique and elegant for those looking for a fashionable option

Fiber Content

- 64% Nylon, 36% Spandex (content can vary slightly by style)
- Latex free



20–30 mmHg
Closed Toe Thigh Medium Long in Graphite

712NMLW91



Sensinnov Scallopped for Onyx Stripes and Sensinnov band on thigh



Scan here to see product images online
sigvaris.com




SIZING CHART

	Small		Medium		Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")	
Thigh Circumference	45.5–66 cm (18"–26")		48–73 cm (19"–28.5")		52–84 cm (20.5"–33")	
Leg Length	Short	Long	Short	Long	Short	Long
	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")
Hip Circumference	83–112.5 cm (32.5"–44.5")		87–127 cm (34.5"–50")		98–150.5 cm (38.5"–59.5")	
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

STYLES AND COLOURS

Closed Toe = ●



Thigh w/Grip-Top (N) Sensinnov grip-top band	●	●	●
Pantyhose (P)	●	●	●

Product Code Example: 712NLSW19

Series	mmHg Level	Style	Size	Gender	Colour
71	2	N	LS	W	19



SEA ISLAND COTTON 220

for women & men

Compression Level

20–30 mmHg

Also available in 15–20mmHg, see page 60

Features and Benefits



Soft



Cotton



Breathable

- Knitted with the world’s rarest cotton for a soft, luxurious feel
- Ideal for sensitive skin
- Breathable and comfortable
- Suitable for those who prefer natural fibers for everyday wear
- Smooth, flat toe seam

Fiber Content

- 43% Sea Island Cotton, 39% Nylon, 11% Spandex, 7% Supima Cotton
- Latex free



Scan here to see product images online
sigvaris.com



20–30 mmHg
Closed Toe Calf Small Long in Brown

222CSLW11



20–30 mmHg
Closed Toe Calf Medium Long in Black

222CMLM99

SIZING CHART




	Small		Medium		Large		X Large	
Ankle Circumference	18–21.5 cm (7"-8.5")		21.5–25.5 cm (8.5"-10")		25.5–29 cm (10"-11.5")		29–36 cm (11.5"-14")	
Calf Circumference	28–39 cm (11"-15.5")		31.5–45 cm (12.5"-17.5")		35.5–51 cm (14"-20")		41–61 cm (16"-24")	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Size	SS	SL	MS	ML	LS	LL	XS	XL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS

Closed Toe = ●



Calf (C)	●	●	●
----------	---	---	---

Product Code Example: 222CLLM10

Series	mmHg Level	Style	Size	Gender	Colour
22	2	C	LL	M	10

MERINO WOOL 240

for women & men

Compression Level

20–30 mmHg

Also available in 15–20 mmHg, see page 62

Features and Benefits



Wool



Thermoregulating



Stylish

- High quality extra-fine Merino wool makes this sock soft and breathable
- Moisture wicking capabilities keep feet dry
- Naturally odour resistant
- Thermoregulating for year-round wear
- Smooth, flat toe seam

Fiber Content

- 53% Extra-fine Merino Wool, 38% Nylon, 9% Spandex
- Latex free



Scan here to see product images online
sigvaris.com

20–30 mmHg
Closed Toe Calf Medium Long in Brown

242CMLM11

20–30 mmHg
Closed Toe Calf Medium Long in Navy

242CMLW10

SIZING CHART


	Small		Medium		Large		X Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")		29–36 cm (11.5"–14")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")		41–61 cm (16"–24")	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Size	SS	SL	MS	ML	LS	LL	XS	XL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12


Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS


Closed Toe = ●



NAVY
10



BROWN
11



BLACK
99

Calf (C)	●	●	●
----------	---	---	---

Product Code Example: 242CMSW99

Series	mmHg Level	Style	Size	Gender	Colour
24	2	C	MS	W	99



LINEN 250

for women & men

Compression Level

20–30 mmHg

Also available in 15–20 mmHg, see page 64

Features and Benefits

Hypoallergenic
Natural fibers

Anti-bacterial

Moisture
Management

Strong and
Durable

- Lightweight and breathable for all day comfort
- Ideal for sensitive skin with hypoallergenic properties
- Strong, stylish and sustainable in four unique colours

Fiber Content

- 48% Cotton, 23% Linen, 22% Nylon and 7% Spandex
- Latex free



20–30 mmHg
Closed Toe Calf Medium Long
in Brown
252CMLM111



20–30 mmHg
Closed Toe Calf Medium Long
in Denim
252CMLW201



Scan here to see product images online
sigvaris.com

SIZING CHART

	Small		Medium		Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")	
Calf Length	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12
Men's Shoe sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS

Closed Toe = ●

BROWN
111

DENIM
201

LAVENDER
207

LIGHT GREY
205

Calf (C)	●	●	●	●
----------	---	---	---	---

Product Code Example: 252CMLW201

Series	mmHg Level	Style	Size	Gender	Colour
25	2	C	ML	W	201



MICROFIBER 820

for men

Compression Level

15–20 mmHg

20–30 mmHg

30–40 mmHg

Features and Benefits



Microfiber



Soft



Stylish

- 100% of microfiber on inside of sock to improve comfort, softness, and facilitate moisture wicking.
- Glides on and off with ease
- Comfortable and durable sock for professional and everyday wear
- Available in extended foot sizes
- Smooth, flat toe seam

Fiber Content

- 87% Nylon , 13% Spandex
(content can vary slightly by style)
- Latex free



30–40 mmHg
Closed Toe Thigh Medium Long
in Black

823NMLM99



Silicone dotted grip-top option
available for calf and designed with
SensinnoV grip-top for thigh



Scan here to see product images online
sigvaris.com

SIZING CHART

	Small		Medium		Large		X Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")		29–36 cm (11.5"–14")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")		41–61 cm (16"–24")	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Thigh Circumference	45.5–66 cm (18"–26")		48–73 cm (19"–28.5")		52–84 cm (20.5"–33")		NA	
Leg Length	Short	Long	Short	Long	Short	Long		
	≤76 cm (30")	>76 cm (30")	≤76 cm (30")	>76 cm (30")	≤76 cm (30")	>76 cm (30")		
Size Shoe Size	SS 7–11	SL 9.5–14	MS 7–11	ML 9.5–14	LS 7–11	LL 9.5–14	XS 7–11	XL 9.5–14
				MX* 12.5–17.5		LX* 12.5–17.5		XX* 12.5–17.5

Men's shoe sizes: Short Leg 7–11, Long Leg 9.5–14 (For extended foot sizes: long leg 12.5–17.5)

*MX, LX, and XX extended foot sizes only available in calf style in Black & Tan Khaki in 15–20 & 20–30 mmHg only. Extended foot sizes not available with grip top

STYLES AND COLOURS

Closed Toe = ●	<div><div>TAN KHAKI 32</div><div>STEEL GREY 93***</div><div>BLACK 99</div></div>		
Calf (C)	●	●	●
Calf w/Grip-Top (C/S)** Silicone dotted grip-top band	●		●
Thigh w/Grip-Top (N) SensinnoV grip-top band	●		●

** Calf w/ Grip-Top not available in 15–20 mmHg

*** Steel Grey only available in 15–20 mmHg and 20–30 mmHg

Product Code Example: 822NLSM99

Series	mmHg Level	Style	Size	Gender	Colour
82	2	N	LS	M	99

Product Code Example: 822CMLM32/S

Series	mmHg Level	Style	Size	Gender	Colour	Grip-top
82	2	C	ML	M	32	/S



MICROFIBER PATTERNS 830

for women & men

Compression Level

20–30 mmHg

Also available in 15–20 mmHg, see page 66

Features and Benefits



Microfiber



Stylish



Patterns

- 100% of microfiber on inside of sock to improve comfort, softness, and facilitate moisture wicking.
- Fashionable, durable sock for everyday wear
- Glides on and off with ease
- Smooth, flat toe seam

Fiber Content

- 87% Nylon, 13% Spandex
(content can vary slightly by style)
- Latex free



Scan here to see product images online
sigvaris.com



20–30 mmHg
Closed Toe Calf Small Long in Royal Blue Argyle

832CSLM48



20–30 mmHg
Closed Toe Calf Large Short in Pink Stripe

832CLSW15

SIZING CHART

	Small		Medium		Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")	
Calf Length	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

Men's Shoe sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS

Closed Toe = ●							
	MARINER STRIPE 02*	GRAPHITE MINI-STRIPE 07	DARK NAVY STRIPE 13	PINK STRIPE 15*	GRAPHITE STRIPE 17	ONYX STRIPE 19	GRAPHITE HEATHER 27
Calf (C)	●	●	●	●	●	●	●
	GREY PLAID 22*	BLACK PLAID 42**	PINK ARGYLE 45*	PURPLE ARGYLE 46	GRAPHITE ARGYLE 47	ROYAL BLUE ARGYLE 48	ONYX ARGYLE 49
Calf (C)	●	●	●	●	●	●	●

* Available for women only

** Available for men only

Product Code Example: 832CMLM17

Series	mmHg Level	Style	Size	Gender	Colour
83	2	C	ML	M	17



Push Yourself

Why choose **Motion**?

- Cushioned sole for added comfort
- Durable thermoregulating wool (422)
- Moisture-wicking yarns

MOTION PRODUCTS

Series	Product	Page
360	Cushioned Cotton	38
412	High Tech	40
422	Thermoregulating Wool	42
850	Comfort	44

BEST SUITED FOR

- Management of Chronic Venous Disease (CVD)
- People who stand for prolonged periods of time
- Casual and occupational footwear
- Walking and everyday activities



CUSHIONED COTTON 360

for women & men

Compression Level

20–30 mmHg

Also available in 15–20 mmHg, see page 58

Features and Benefits



Cushioned Sole



Cotton



Sensitive Skin

- Supima cotton for highest level of comfort and breathability
- Extra-padded sole, heel, and toe to protect feet
- Ideal for sensitive skin and allergies
- Flat, low friction toe seam
- Diabetic friendly
- Ideal for people who stand for prolonged periods of time

Fiber Content

- 64% Supima Cotton, 28% Nylon, 8% Spandex
- Latex free



Scan here to see product images online
sigvaris.com



20–30 mmHg
Closed Toe Calf Medium Long in Black

362CMLW99



Flat, low friction toe seam


SIZING CHART

	Small		Medium		Large		X Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")		29–36 cm (11.5"–14")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")		41–61 cm (16"–24")	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Size	SS	SL	MS	ML	LS	LL	XS	XL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12
Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS

Closed Toe = ●



Calf (C)	●	●
----------	---	---

Product Code Example: 362CMLM00

Series	mmHg Level	Style	Size	Gender	Colour
36	2	C	ML	M	00



HIGH TECH 412

for women & men

Compression Level

20–30 mmHg

Features and Benefits



Moisture
Wicking



Anatomically
Contoured



Colours

- High-Tech fibers for moisture wicking and improved thermoregulation
- Cushioned footbed, toe, and Achilles tendon zones for optimal comfort
- Anatomically contoured to support the natural shape of the foot and calf
- Smooth, flat toe seam
- Perfect for everyday activities and active individuals

Fiber Content

- 64% Nylon, 17% Polyester, 19% Spandex
- Latex free




Scan here to see product images online
sigvaris.com



20–30 mmHg
Closed Toe Calf Small Medium in
Limeade

412CSM64



20–30 mmHg
Closed Toe Calf Medium Large in
Pink

412CML56

SIZING CHART

	Small			Medium			Large			X Large		
Ankle Circumference	18–21.5 cm (7"–8.5")			21.5–25.5 cm (8.5"–10")			25.5–29 cm (10"–11.5")			29–33 cm (11.5"–13")		
Calf Circumference	28–39 cm (11"–15.5")			31.5–45 cm (12.5"–17.5")			35.5–51 cm (14"–20")			40.5–58 cm (16"–23")		
Women's Shoe Size	4.5–7.5	8–10.5	11–13.5	4.5–7.5	8–10.5	11–13.5	8–10.5	11–13.5	>14	8–10.5	11–13.5	>14
Men's Shoe Size	≤6.5	7–9.5	10–12.5	≤6.5	7–9.5	10–12.5	7–9.5	10–12.5	13–15.5	7–9.5	10–12.5	13–15.5
Size	SS	SM	SL	MS	MM	ML	LS	LM	LL	XS	XM	XL

STYLES AND COLOURS

Closed Toe = ●	WHITE 00	BLUE 50	LIME 54	PINK 56	RED 59	STEEL BLUE 62	BLOOD ORANGE 63	LIMEADE 64	BLACK 99
	●	●	●	●	●	●	●	●	●

Product Code Example: 412CML62

Series	mmHg Level	Style	Size	Colour
41	2	C	ML	62



THERMO-REGULATING WOOL 422

for women & men

Compression Level

20–30 mmHg

Also available in 15–20 mmHg, see page 76

Features and Benefits



Wool



Cushioned Sole



Thermoregulating

- Natural properties of wool make these socks odour resistant
- Thermoregulating and moisture-wicking for year-round wear
- Cushioned, padded soles help prevent blisters
- Soft and comfortable
- Smooth, flat toe seam
- Perfect for outdoor work environments and activities

Fiber Content

- 67% Fine Merino Wool, 26% Nylon, 7% Spandex
- Latex free



Scan here to see product images online
sigvaris.com

20–30 mmHg
Closed Toe Calf Large Long
in Charcoal

422CLL12

20–30 mmHg
Closed Toe Calf Medium Short
in Olive


422CMS31

SIZING CHART

	Small			Medium			Large		
Ankle Circumference	18–21.5 cm (7"–8.5")			21.5–25.5 cm (8.5"–10")			25.5–29 cm (10"–11.5")		
Calf Circumference	28–39 cm (11"–15.5")			31.5–45 cm (12.5"–17.5")			35.5–51 cm (14"–20")		
Women's Shoe Size	4.5–7.5	8–10.5	11–13.5	4.5–7.5	8–10.5	11–13.5	8–10.5	11–13.5	>14
Men's Shoe Size	≤6.5	7–9.5	10–12.5	≤6.5	7–9.5	10–12.5	7–9.5	10–12.5	13–15.5
Size	SS	SM	SL	MS	MM	ML	LS	LM	LL

STYLES AND COLOURS

Closed Toe = ●



Calf (C)	●	●
----------	---	---

Product Code Example: 422CML12

Series	mmHg Level	Style	Size	Colour
42	2	C	ML	12



COMFORT 850

for women & men

Compression Level

20–30 mmHg

Features and Benefits



Cushioned Sole



Durable



Moisture Wicking

- NILIT Aquarius fibers provide moisture wicking to keep feet dry and comfortable
- Cushioned soles to absorb impact and help prevent blisters
- Smooth, flat toe seam
- Suitable for uniforms requiring solid black color

Fiber Content

- 86% Nylon, 14% Spandex
- Latex free



20–30 mmHg
Closed Toe Calf Large Long in Black

852CLLM99



20–30 mmHg
Closed Toe Calf Medium Long in Black

852CMLW99



Scan here to see product images online
sigvaris.com

SIZING CHART


	Small		Medium		Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")	
Calf Length	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

Men's Shoe sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS

Closed Toe = ●



Calf (C)	●
----------	---

Product Code Example: 852CMLM99

Series	mmHg Level	Style	Size	Gender	Colour
85	2	C	ML	M	99



Take care of Yourself

Why choose **Specialty**?

- Specific medical needs
- Need low fatigue fibers
- Natural rubber (500)

SPECIALTY PRODUCTS

Series	Product	Page
500	Natural Rubber	48
910	Advance Armsleeve	50

BEST SUITED FOR

- Severe varicosities
- Chronic edema post DVT
- Severe venous stasis
- Early stages of lymphedema
- Active venous ulcers and prevention of ulcer recurrence
- People with early stages of arm lymphedema
- Post-mastectomy care



NATURAL RUBBER 500

for women & men

Compression Level

30–40 mmHg

40–50 mmHg

50–60 mmHg

Features and Benefits



Natural Rubber



Durable



Largest Size Range

- Durable, long-term containment solution with minimal fiber fatigue throughout the day
- Extended sizing for a custom-like fit
- **Ideal for patients with more serious indications, including:**
 - Severe varicosities
 - Chronic edema post DVT
 - Prevention and treatment of Active Venous Ulcers
 - Severe venous stasis
 - Lymphedema

Fiber Content

- 70% Nylon, 30% Natural Latex Rubber (Calf & Thigh-high styles)
- 75% Nylon, 20% Natural Latex Rubber, 5% other fibers (Pantyhose style)
- **Caution:** contains natural rubber latex that may cause allergic reaction



30–40 mmHg
Open Toe Thigh Small Average Long
in Beige

503NS2O77



Silicone dotted grip-top available
on thigh

SIZING CHART

	Small				Medium				Large				X Large			
Ankle Circumference	19–23 cm (7.5"–9")				23–30 cm (9"–12")				30–40 cm (12"–15.5")				40–46 cm (15.5"–18")			
Calf Circumference	Average		Full		Average		Full		Average		Full		Average		Full	
	28–38 cm (11"–15")		≤43 cm (17")		32–43 cm (12.5"–17")		≤48 cm (19")		36–48 cm (14.5"–19")		≤57 cm (22.5")		45–57 cm (18"–22.5")		≤68 cm (27")	
Calf Circumference for 505C (Full only)	505C: full calf only		32–43 cm (12.5"–17")		505C: full calf only		36–48 cm (14.5"–19")		505C: full calf only		45–57 cm (18"–22.5")		505C: not available in XL			
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long
	≤40 cm (≤15.5")	>40 cm (>15.5")	≤40 cm (≤15.5")	>40 cm (>15.5")	≤40 cm (≤15.5")	>40 cm (>15.5")	≤40 cm (≤15.5")	>40 cm (>15.5")	≤40 cm (≤15.5")	>40 cm (>15.5")	≤40 cm (≤15.5")	>40 cm (>15.5")	≤40 cm (≤15.5")	>40 cm (>15.5")	≤40 cm (≤15.5")	>40 cm (>15.5")
Thigh Circumference	Average		Full		Average		Full		Average		Full		N/A			
	41–54 cm (16"–21")		≤61 cm (24")		48–61 cm (19"–24")		≤71 cm (28")		58–71 cm (23"–28")		≤81 cm (32")					
Thigh Circumference for 505T (Full only)	505T: full thigh only		≤61 cm (24")		505T: full thigh only		≤71 cm (28")		505T: full thigh only		≤81 cm (32")					
Leg Length	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long				
	≤71 cm (28")	>71 cm (28")	≤71 cm (28")	>71 cm (28")	≤71 cm (28")	>71 cm (28")	≤71 cm (28")	>71 cm (28")	≤71 cm (28")	>71 cm (28")	≤71 cm (28")	>71 cm (28")				
Hip Circumference	<122 cm (48")				<137 cm (54")				<155 cm (61")							
Size	S1	S2	S3	S4	M1	M2	M3	M4	L1	L2	L3	L4	X1	X2	X3	X4

STYLES AND COLOURS

Open Toe = ○



Caution: this product contains natural rubber latex which may cause allergic reactions. For all indications, cautions, and contraindications, please see page 101 for details. Sigvaris accepts no liability for non-observance of contraindications or cautions.

Calf (C)	○	Calf only available in S, M, and L full circumference for 50–60 mmHg
Thigh w/out Grip-Top (T)	○	
Thigh w/Grip-Top (N) Silicone dotted grip-top band	○	Not available in 50-60 mmHg or XL
Pantyhose (P)	○	Pantyhose only available in 30–40 mmHg and only in full thigh circumference
Thigh w/Waist Attachment (W) Indicate Leg LEFT (/L) or RIGHT (/R)	○	Thigh w/ waist attachment is available in average and full thigh circumferences in 30–40 mmHg and in full thigh circumferences only in 40–50 mmHg

Product Code Example: 504CS2O77

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour
50	4	C	S2	O	77

Product Code Example: 503WL1O77/L

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour	Leg
50	3	W	L1	O	77	/L



ADVANCE ARMSLEEVE 910

for women & men

Compression Level

- 15–20 mmHg
- 20–30 mmHg
- 30–40 mmHg

Features and Benefits



Moisture Wicking



Sensinnov Grip-top Band



Breathable

- Ideal for patients with arm lymphedema or post-mastectomy care
- Provides moisture wicking properties to keep skin dry and comfortable
- Breathable and soft
- Sensinnov grip-top band ensures armsleeve stays in place all day long
- Available with or without gauntlet

Fiber Content

- 50% Nylon, 28% Rayon, 22% Spandex
- Latex free



With Gauntlet

20–30 mmHg
Armsleeve With Gauntlet
Small Long in Beige
912ASLG77/S



Sensinnov grip-top band

SIZING CHART

	Small		Medium				Large				X Large			
Palm (A)	17–21 cm (6.5"–8")		18–22 cm (7"–8.5")				19–23 cm (7.5"–9")				20–24 cm (8"–9.5")			
Wrist (C)	14–16 cm (5.5"–6.5")		16–18 cm (6.5"–7")				18–20 cm (7"–8")				20–23 cm (8"–9")			
Upper Mid-Arm (F)	Regular		Regular		Plus		Regular		Plus		Regular		Plus	
	26–34 cm (10"–13.5")		29–37 cm (11.5"–14.5")		37–43 cm (14.5"–17")		32–40 cm (12.5"–15.5")		40–46 cm (15.5"–18")		34–42 cm (13.5"–16.5")		42–48 cm (16.5"–19")	
Arm Length (C-G)*	Regular	Long	Regular	Long	Regular	Long	Regular	Long	Regular	Long	Regular	Long	Regular	Long
	≤43 cm (17")	>43 cm (17")	≤43 cm (17")	>43 cm (17")	≤43 cm (17")	>43 cm (17")	≤43 cm (17")	>43 cm (17")	≤43 cm (17")	>43 cm (17")	≤43 cm (17")	>43 cm (17")	≤43 cm (17")	>43 cm (17")
Size	SR	SL	MR	ML	MR/P	ML/P	LR	LL	LR/P	LL/P	XR	XL	XR/P	XL/P


*Arm should have slight bend in arm when measuring length

STYLES AND COLOURS

Open Sleeve = ●

	BEIGE 77
Armsleeve w/Gauntlet (G) Sensinnov grip-top band	●
Armsleeve w/out Gauntlet (O) Sensinnov grip-top band	●

Note: all armsleeves come with Grip-top (/S)
15–20 mmHg not available with gauntlet or plus sizing



ARMSLEEVE WITHOUT GAUNTLET

Product Code Example: 912ASLO77/S

Series	mmHg Level	Style	Size	W/out Gauntlet	Colour	Grip-Top
91	2	A	SL	O	77	/S

ARMSLEEVE WITH GAUNTLET

Product Code Example: 912ALRG77/PS

Series	mmHg Level	Style	Size	W/ Gauntlet	Colour	Plus-Size	Grip-Top
91	2	A	LR	G	77	/P	S



MARINER STRIPE

Care for Yourself

Why choose **Consumer Health**?

- Tired, heavy, restless, aching legs
- Travelers, occupational leg health risks, expecting mothers, family history of Chronic Venous Disease, prevention of swelling and varicose veins
- Simplified sizing

CONSUMER HEALTH PRODUCTS

Series	Product	Page
146 / 186	Casual Cotton	54
146/T / 186/T	Traveno Travel Socks	56
142 / 182	Cushioned Cotton	58
151 / 191	Sea Island Cotton	60
152 / 192	All-Season Merino Wool	62
153 / 193	Linen	64
143 / 183	Microfiber Shades	66
120	Sheer Fashion	68
170	Soft Silhouette Leggings	70
401	Athletic Recovery Socks	72
412	Performance Sleeves	74
421	Merino Outdoor Socks	76
160	Eversoft Diabetic Socks	78



CASUAL COTTON

146/186

for women & men

Compression Level

15–20 mmHg

Features and Benefits

- Made with premium Supima Cotton for optimum softness and durability
- Subtle ribbed design
- Flat, low-friction toe seam to prevent irritation

Fiber Content

- 66% Supima Cotton, 26% Nylon, 8% Spandex
- Latex free

SIZING CHART				
146	Women's Shoe	5–7	7.5–9.5	10–12
	Ankle Circumference	18–23 cm (7"–9")	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")
186	Men's Shoe	6.5–8.5	9–11	11.5–14
	Ankle Circumference	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")	21–29 cm (8.5"–11.5")
Size		A	B	C

STYLES AND COLOURS					
Closed Toe = ●					
	WHITE 00	NAVY 10	BROWN 11	KHAKI 30	BLACK 99
Calf (C)	●	●	●	●	●

Product Code Example: 146CB10


Series	Style	Size	Colour
146	C	B	10



146
Women's Calf
Style

15–20 mmHg
Closed Toe Calf Size B in Khaki

146CB30



186
Men's Calf
Style

15–20 mmHg
Closed Toe Calf Size C in Black

186CC99



TRAVENO TRAVEL SOCKS 146/T - 186/T

for women & men

Compression Level

15–20 mmHg

Features and Benefits

- Made with premium Supima Cotton for optimum softness and durability
- Subtle ribbed design
- Flat, low-friction toe seam to prevent irritation

Fiber Content

- 66% Supima Cotton, 26% Nylon, 8% Spandex
- Latex free

STYLES AND COLOURS

Closed Toe = ●



BLACK
99

Calf (C)



SIZING CHART

146/T	Women's Shoe	5–7	7.5–9.5	10–12
	Ankle Circumference	18–23 cm (7"–9")	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")
186/T	Men's Shoe	6.5–8.5	9–11	11.5–14
	Ankle Circumference	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")	21–29 cm (8.5"–11.5")
	Size	A	B	C

IMPROVE YOUR
CIRCULATION
IN THE AIR
OR ON
THE GROUND

Traveno
TRAVEL
SOCKS



Travel Socks for Women & Men
Available in 15–20 mmHg in black

Countertop displays available. Contact your Territory Manager or call **1.800.363.4999** for more information

CUSHIONED COTTON 142/182

for women & men

Compression Level

15–20 mmHg

Also available in 20–30 mmHg, see page 38

Features and Benefits

- Supima cotton for highest level of comfort and breathability
- Extra-padded sole, heel, and toe to protect feet
- Ideal for sensitive skin and allergies
- Flat, low friction toe seam
- Diabetic friendly
- Ideal for people who stand for prolonged periods of time

Fiber Content

- 68% Supima Cotton, 25% Nylon, 7% Spandex
- Latex free

SIZING CHART

142	Women's Shoe	5–7	7.5–9.5	10–12
	Ankle Circumference	18–23 cm (7"–9")	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")
182	Men's Shoe	6.5–8.5	9–11	11.5–14
	Ankle Circumference	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")	21–29 cm (8.5"–11.5")
	Size	A	B	C

STYLES AND COLOURS

Closed Toe = ●

	WHITE 00	BLACK 99
Calf (C)	●	●

Product Code Example: 142CB00

Series	Style	Size	Colour
142	C	B	00



142
Women's Calf Style

15–20 mmHg
Closed Toe Calf Size B in White

142CB00



182
Men's Calf Style

15–20 mmHg
Closed Toe Calf Size C in Black

182CC99



SEA ISLAND COTTON 151/191

for women & men

Compression Level

15–20 mmHg

Also available in 20–30 mmHg, see page 26

Features and Benefits

- Knitted with the world’s rarest cotton for a soft, luxurious feel
- Ideal for sensitive skin
- Breathable and comfortable
- Suitable for those who prefer natural fibers for everyway wear
- Smooth, flat toe seam

Fiber Content

- 43% Sea Island Cotton, 39% Nylon, 11% Spandex, 7% Supima Cotton
- Latex free

SIZING CHART				
151	Women's Shoe	5–7	7.5–9.5	10–12
	Ankle Circumference	18–23 cm (7"–9")	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")
191	Men's Shoe	6.5–8.5	9–11	11.5–14
	Ankle Circumference	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")	21–29 cm (8.5"–11.5")
	Size	A	B	C

STYLES AND COLOURS			
Closed Toe = ●	<div><div>NAVY 10</div><div>BROWN 11</div><div>BLACK 99</div></div>		
	Calf (C)	●	●

Product Code Example: 151CA10

Series	Style	Size	Colour
151	C	A	10



15–20 mmHg
Closed Toe Calf Size A in Navy

151CA10



15–20 mmHg
Closed Toe Calf Size B in Brown

191CB11



ALL-SEASON MERINO WOOL 152/192

for women & men

Compression Level

15–20 mmHg

Also available in 20–30 mmHg, see page 28

Features and Benefits

- High quality extra-fine Merino Wool makes this sock soft and breathable
- Moisture wicking capabilities and naturally odour resistant
- Thermoregulating for year-round wear and keeps feet dry
- Smooth, flat toe seam

Fiber Content

- 53% Extra-fine Merino Wool, 38% Nylon, 9% Spandex
- Latex free

SIZING CHART				
152	Women's Shoe	5–7	7.5–9.5	10–12
	Ankle Circumference	18–23 cm (7"–9")	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")
192	Men's Shoe	6.5–8.5	9–11	11.5–14
	Ankle Circumference	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")	21–29 cm (8.5"–11.5")
	Size	A	B	C

STYLES AND COLOURS			
Closed Toe = ●	<div><div>NAVY 10</div><div>BROWN 11</div><div>BLACK 99</div></div>		
	Calf (C)	●	●

Product Code Example: 192CB11

Series	Style	Size	Colour
192	C	B	11



15–20 mmHg
Closed Toe Calf Size A in Black

152CA99



15–20 mmHg
Closed Toe Calf Size C in Navy

192CC10



LINEN 153/193

for women & men

Compression Level

15–20 mmHg

Also available in 20–30 mmHg, see page 30

Features and Benefits

- Lightweight and breathable for all day comfort
- Ideal for sensitive skin with hypoallergenic properties
- Strong, stylish and sustainable in four unique colours

Fiber Content

- 48% Cotton, 23% Linen, 22% Nylon and 7% Spandex
- Latex free

SIZING CHART				
153	Women's Shoe	5–7	7.5–9.5	10–12
	Ankle Circumference	18–23 cm (7"–9")	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")
193	Men's Shoe	6.5–8.5	9–11	11.5–14
	Ankle Circumference	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")	21–29 cm (8.5"–11.5")
	Size	A	B	C

STYLES AND COLOURS				
Closed Toe = ●				
		BROWN 111	DENIM 201	LAVENDER 207
Calf (C)		●	●	●
		LIGHT GREY 205		

Product Code Example: 153CA207

Series	Style	Size	Colour
153	C	A	207



153
Women's
Calf Style

15–20 mmHg
Closed Toe Calf Size A in Brown

153CA111



193
Men's
Calf Style

15–20 mmHg
Closed Toe Calf Size C in Denim

193CC201



MICROFIBER SHADES 143/183

for women & men

Compression Level

15–20 mmHg

Also available in 20–30 mmHg, see page 34

Features and Benefits

- 100% of microfiber sits next to skin for a soft feel
- Fashionable, durable socks for everyday wear
- Excellent breathability and moisture management
- Glides on and off with ease
- Smooth, flat toe seam

Fiber Content

- 87% Nylon, 13% Spandex
- Latex free



143
Women's
Calf Style

15–20 mmHg
Closed Toe Calf Size B in Pink Argyle

143CB45



183
Men's
Calf Style

15–20 mmHg
Closed Toe Calf Size C in Black Plaid

183CC42

SIZING CHART

143	Women's Shoe	5–7	7.5–9.5	10–12
	Ankle Circumference	18–23 cm (7"–9")	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")
183	Men's Shoe	6.5–8.5	9–11	11.5–14
	Ankle Circumference	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")	21–29 cm (8.5"–11.5")
	Size	A	B	C

STYLES AND COLOURS

Closed Toe = ●

Actual colour may vary

Calf (C)	●	●	●	●	●	●	●
Calf (C)	●	●	●	●	●	●	●

* Available for women only

** Available for men only

Product Code Example: 143CB02

Series	Style	Size	Colour
143	C	B	02



SHEER FASHION 120

for women

Compression Level

15–20 mmHg

Features and Benefits

- Incredible sheer design in a variety of fashionable colours
- Sheer toe and heel area
- Comfortable stretch and soft fabric, excellent for everyday wear

Fiber Content

- 79% Nylon, 21% Spandex
(Content can vary slightly by style)
- Latex free



15–20 mmHg
Closed Toe Calf Size C in Charcoal

120CC12



15–20 mmHg
Closed Toe Pantyhose Size C in Black

120PC99

PANTYHOSE SIZING CHART

cm/ft	Kg/Lbs	45 / 100	48 / 105	50 / 110	52 / 115	54 / 120	57 / 125	59 / 130	61 / 135	64 / 140	66 / 145	68 / 150	70 / 155	73 / 160	75 / 165	77 / 170	80 / 175	82 / 180	84 / 185	86 / 190	89 / 195	91 / 200	93 / 205	95 / 210	98 / 215	100 / 220	102 / 225	105 / 230
145 cm / 4'9"																												
147 cm / 4'10"																												
150 cm / 4'11"																												
152 cm / 5'0"																												
155 cm / 5'1"																												
157 cm / 5'2"																												
160 cm / 5'3"																												
162 cm / 5'4"																												
165 cm / 5'5"																												
167 cm / 5'6"																												
170 cm / 5'7"																												
173 cm / 5'8"																												
175 cm / 5'9"																												
178 cm / 5'10"																												
180 cm / 5'11"																												
183 cm / 6'0"																												

CALF / THIGH SIZING CHART

120	Women's Shoe	5–7	7.5–9.5	10–12
	Ankle Circumference	18–23 cm (7"–9")	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")
	Size	A	B	C

Weight Guidelines For Thigh Style: Size A: <135Lbs ; Size B: <180Lbs ; Size C: <230Lbs.

STYLES AND COLOURS

Closed Toe = ●

Open Toe = ○

	CHARCOAL 12	TAUPE 29	HONEY 33	GOLDEN 36	CAFÉ 73	BLACK 99
Calf (C)	●	●○	●○	●○	●	●
Thigh (N)	●	●	●	●	●	●
Pantyhose (P)	●	●	●	●	●	●

Product Code Example: 120CCO36

Series	Style	Size	Open Toe	Colour
120	C	C	O	36

Product Code Example: 120PE73

Series	Style	Size	Colour
120	P	E	73



SOFT SILHOUETTE LEGGINGS 170

for women

Compression Level

15–20 mmHg

Features and Benefits

- Advanced soft stretch microfiber fabric shapes legs and smooths skin
- Foldable, comfort-stretch waist band provides a seamless look for all-day wearing comfort
- Stylish, subtle ribbed design looks chic and pairs perfectly with dresses, blouses, tunics, and shirts

Fiber Content

- 83% Nylon, 17% Spandex
- Latex free



15–20 mmHg
Legging Size B in Black

170LB99



Non-binding ankle cuff
stays comfortably in place
for all day wear

SIZING CHART

cm / ft	Kg / Lbs	45 / 100	48 / 105	50 / 110	52 / 115	54 / 120	57 / 125	59 / 130	61 / 135	64 / 140	66 / 145	68 / 150	70 / 155	73 / 160	75 / 165	77 / 170	80 / 175	82 / 180	84 / 185	86 / 190	89 / 195	91 / 200	93 / 205	95 / 210	98 / 215	100 / 220	102 / 225	105 / 230
147 cm / 4'10"																												
150 cm / 4'11"																												
152 cm / 5'0"																												
155 cm / 5'1"																												
157 cm / 5'2"																												
160 cm / 5'3"																												
162 cm / 5'4"																												
165 cm / 5'5"																												
167 cm / 5'6"																												
170 cm / 5'7"																												
173 cm / 5'8"																												
175 cm / 5'9"																												
178 cm / 5'10"																												
180 cm / 5'11"																												
183 cm / 6'0"																												

Guidelines for Ankle Circumference
Size A: 18–23 cm (7"–9"), Size B & C: 20–26 cm (8"–10"), Size D: 22–29 cm (8.5"–11.5")

STYLES AND COLOURS

Open Foot = 	 BLACK 99
Leggings (L)	

Product Code Example: 170LB99

Series	Style	Size	Colour
170	L	B	99



ATHLETIC
RECOVERY
SOCKS 401

for women & men

Compression Level

15–20 mmHg

Features and Benefits

- Drirelease yarns provide moisture wicking and odour control
- Comfortable and easy to put on and take off
- Helps reduce lactic acid buildup
- Mild cushioning in sole, heel, and toe area for added comfort

Fiber Content


- 67% drirelease (57% Polyester, 10% Cotton), 26% Nylon, 7% Spandex
- Latex free

SIZING CHART

401	Women's Shoe	<7.5	7.5–9.5	10–12	12.5–15	≥15.5
	Men's Shoe	<6.5	6.5–8.5	9–11	11.5–14	14.5–17
	Ankle Circumference	18–24 cm (7"–9.5")	19–26 cm (7.5"–10")	20–27 cm (8"–10.5")	21–29 cm (8.5"–11.5")	21–29 cm (8.5"–11.5")
	Size	S	M	L	X	2

STYLES AND COLOURS

Closed Toe = ●



Calf (C)	●	●
----------	---	---

Product Code Example: 401C200

Series	Style	Size	Colour
401	C	2	00



15–20 mmHg
Closed Toe Calf Size X-Large in Black

401CX99



15–20 mmHg
Closed Toe Calf Size Medium in White

401CM00



PERFORMANCE SLEEVES 412

for women & men

Compression Level

20–30 mmHg

Features and Benefits

- Stabilizes and supports active muscles
- Ideal for running, triathlons, endurance sports, and gym workouts
- Fabric resists deterioration in salt, chlorine, or fresh water
- Quick drying
- High-visibility reflective transfer

Fiber Content

- 75% Nylon, 25% Spandex
- Latex free

SIZING CHART					
412	Ankle Circumference	18–21.5 cm (7"–8.5")	21.5–25.5 cm (8.5"–10")	25.5–29 cm (10"–11.5")	29–33 cm (11.5"–13")
	Calf Circumference	28–39 cm (11"–15.5")	31.5–45 cm (12.5"–17.5")	35.5–51 cm (14"–20")	40.5–58 cm (16"–23")
Size		S	M	L	X

STYLES AND COLOURS						
		<div></div> <div>WHITE 00</div>	<div></div> <div>BLUE 50</div>	<div></div> <div>PINK 56*</div>	<div></div> <div>PURPLE 58*</div>	<div></div> <div>BLACK 99</div>
Calf Sleeve = <input type="radio"/>						
Calf Sleeve (V)		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

* Pink and purple not available in X-Large

Product Code Example: 412VL56

Series	Style	Size	Colour
412	V	L	56



20–30 mmHg
Calf Sleeve Size Small in Purple

412VS58



20–30 mmHg
Calf Sleeve Size Medium in Blue

412VM50

MERINO OUTDOOR SOCKS 421

for women & men

Compression Level

15–20 mmHg

Also available in 20–30 mmHg, see page 42

Features and Benefits

- Natural properties of wool make these socks odour resistant
- Thermoregulating and moisture-wicking for year-round wear
- Cushioned, padded soles help prevent blisters
- Mesh flex zones for comfort and ventilation
- Smooth, flat toe seam
- Soft and comfortable
- Perfect for outdoor work environments and activities

Fiber Content

- 67% Fine Merino Wool, 26% Nylon, 7% Spandex
- Latex free

SIZING CHART					
421	Women's Shoe	<7.5	7.5–9.5	10–12	12.5–15
	Men's Shoe	<6.5	6.5–8.5	9–11	11.5–14
	Ankle Circumference	18–24 cm (7"–9.5")	19–26 cm (7.5"–10")	20–27 cm (8"–10.5")	21–29 cm (8.5"–11.5")
	Size	S	M	L	X

STYLES AND COLOURS		
Closed Toe = ●	CHARCOAL 12	OLIVE 31
	Calf (C)	● ●

Product Code Example: 421CS31

Series	Style	Size	Colour
421	C	S	31



15–20 mmHg
Closed Toe Calf Size Large in Olive

421CL31

15–20 mmHg
Closed Toe Calf Size Medium in Charcoal

421CM12



EVERSOFT DIABETIC SOCKS 160

for women & men

Compression Level

8-15 mmHg

Features and Benefits

- Padded sole, heel, and toe for comfort and protection
- Flat, low friction toe seam
- Compression-free foot
- Drirelease yarns provide moisture wicking and odour control
- Suitable for individuals with limited dexterity or lower tolerance for compression
- Ideal for individuals with sensitive skin

Fiber Content

- 67% drirelease (57% Polyester, 10% Cotton), 26% Nylon, 7% Spandex
- Latex free



8-15 mmHg
Closed Toe Calf Size Small in White

160CS00



8-15 mmHg
Closed Toe Calf Size Large in White

160CL00

SIZING CHART

160	Women's Shoe	<7.5	7.5-9.5	10-12	12.5-15
	Men's Shoe	<6.5	6.5-8.5	9-11	11.5-14
	Ankle Circumference	18-24 cm (7"-9.5")	19-26 cm (7.5"-10")	20-27 cm (8"-10.5")	21-29 cm (8.5"-11.5")
	Size	S	M	L	X

STYLES AND COLOURS

Closed Toe = ●



Calf (C)	●
----------	---

Product Code Example: 160CS00

Series	Style	Size	Colour
160	C	S	00



Why Choose Inelastic Wraps?

- Adjustable pressure
- Large variety of solutions

MEDICALLY COMPLEX EDEMA PRODUCTS

Product	Page
Compreflex Wrap Overview	82
Compreflex Standard Calf & Foot	84
Compreflex Liner	86
Compребoot	87
Compreflex Strap Extender	88
Accutab	88
Compreflex Thigh	89
Compreflex Knee	90
Compreknee	91
Coolflex	92
Compreshort	96

The SIGVARIS GROUP line of products for medically complex edema (MCE) are for the treatment and management of conditions impacting the lymphatic and venous circulatory systems. The product portfolio contains ready-to-wear garments for the active treatment and maintenance of lymphedema and for patients managing chronic venous disorders in the lower extremities. Sigvaris is looking towards a brighter future with new therapies for complex venous and lymphatic diseases. Our inelastic wraps deliver maximum therapeutic results with increased patient compliance.

COMPREFLEX

20-50 mmHg



Wrap Details

ACCURATE, MEASURED THERAPEUTIC COMPRESSION FOR SELF-DIRECTED TREATMENT AND SELF-MANAGED CARE.

The Compreflex range of products are low profile compression devices designed to deliver moderate to firm compression to the upper and lower extremities. Each Compreflex device is made of a soft conforming inelastic material, offering low to minimal stretch for proper fit and the gradient application of compression. The Compreflex garments are designed with key features that support the patient, including our patented Accutab® Devices which ensure the accurate application of the prescribed dosage of therapeutic compression, an interior stretch panel for easy donning, a contoured calf for limb shaping, and opposing straps that are rounded for improved fit and patient comfort.

INTELLIGENT DESIGN

Garment features a contoured calf, providing a superior fit

PREMIUM MATERIALS

Soft, conforming material provides a comfortable fit

EASY APPLICATION

With a front stretch panel to secure the garment in place and straps that roll back, the garment is easy to don and adjust, improving compliance

SUPERIOR QUALITY

Each garment is double stitched to ensure durability and is backed by a six-month warranty

Features & Benefits


- High quality inelastic wrap composed of durable breathable inelastic material with moisture wicking capacities
- Contoured design helps the wrap conform to the natural shape of the leg
- Front stretch panel secures the garment in place, while the Velcro straps assist for easy application
- Ability to customize with different liner and boot options

Conditions

- Chronic Venous Disease (CVD) including, but not limited to, ulcer due to venous stasis
- Lymphedema - Active and Maintenance
- Non-specific or post surgical edema

How it works

PROVEN ACCURACY FOR RESTING PRESSURE




Accutab has been rigorously tested according to established guidelines in the International Compression Club Consensus documents. Testing was conducted using PicoPress, which is recognized for its high accuracy and reproducibility

Effective working pressure

- Inelastic compression reduces leg edema and improves venous blood flow
- Inelastic compression can increase arterial perfusion in patients with mild to moderate peripheral arterial disease

Resting Pressure

Donning Compreflex



Resting Pressure

Laying Down Sitting Still

Working pressure

Walking, stair climbing, sit to stand, toe raises or other therapeutic exercises to actuate the calf muscle

Safe, precise resting pressure

- Under the supervision of a healthcare provider, inelastic compression can be used safely for patients with mild to moderate peripheral arterial disease.
- Accutab delivers precise dosages for controlled initial resting pressure when donning a Sigvaris Compreflex garment



Standard CALF & FOOT

20–50 mmHg

For those who require a calf and foot garment Includes:

- Compreflex Wrap
- Non-compressive Basic Liner
- Non attached foot piece

PRODUCT CODE					
Size Regular	Black	Beige	Size Tall	Black	Beige
Small	1401-BKR	1411-BKR	Small	1401-BKT	1411-BKT
Medium	1402-BKR	1412-BKR	Medium	1402-BKT	1412-BKT
Large	1403-BKR	1413-BKR	Large	1403-BKT	1413-BKT
X Large	1404-BKR	1414-BKR	X Large	1404-BKT	1414-BKT
XX Large	1405-BKR	1415-BKR	XX Large	1405-BKT	1415-BKT




Note: • 1401-BKR, 1411-BKR, 1401-BKT, and 1411-BKT includes a Small Regular Boot
• 1402-BKR, 1412-BKR, 1402-BKT, 1412-BKT, 1403-BKR, 1413-BKR, 1403-BKT, and 1413-BKT includes a Medium /Large Regular Boot
• 1404-BKR, 1414-BKR, 1404-BKT, 1414-BKT, 1405-BKR, 1415-BKR, 1405-BKT, and 1415-BKT includes a X-Large / XX-Large Regular Boot

Sizing Chart

COMPREFLEX

CALF SIZING CHART (A–C)					
	Small	Medium	Large	X Large	XX Large
C	29–39 cm	34–44 cm	40–50 cm	46–56 cm	54–64 cm
A	19–26 cm	21–30 cm	26–36 cm	31–41 cm	36–46 cm



G Length
Regular: 30–33.9cm
Tall: 34–42cm

Foot Size:
I: Measure circumference of ankle bend and heel
J: Measure circumference across metatarsal heads
K: Measure **medial** length from heel to 1st metatarsal head

Compreflex Calf options



Standard CALF

20–50 mmHg

For those who require an alternative foot garment or do not require a foot piece

Includes:

- Compreflex Wrap
- Non-compressive Basic Liner

Can be matched with alternative foot pieces (see page 87)



Transition CALF

20–50 mmHg

For those who require a compressive liner for the foot and ankle

Includes:

- Compreflex Wrap
- Transition Liner, which provides 10–15 mmHg in the foot and ankle

Ability to wear regular shoe over liner



Complete CALF

20–50 mmHg

For those who require a compressive liner with compression in the foot and ankle, and graduated compression in the calf

Includes:

- Compreflex Wrap
- Complete Liner, which provides 10-15mmHg in the foot and ankle, and graduated compression in the calf

Ability to wear regular shoe over liner

PRODUCT CODE		
Size Regular	Black	Beige
Small	1401-NF-BKR	1411-NF-BKR
Medium	1402-NF-BKR	1412-NF-BKR
Large	1403-NF-BKR	1413-NF-BKR
X Large	1404-NF-BKR	1414-NF-BKR
XX Large	1405-NF-BKR	1415-NF-BKR
Size Tall		
Small	1401-NF-BKT	1411-NF-BKT
Medium	1402-NF-BKT	1412-NF-BKT
Large	1403-NF-BKT	1413-NF-BKT
X Large	1404-NF-BKT	1414-NF-BKT
XX Large	1405-NF-BKT	1415-NF-BKT

PRODUCT CODE		
Size Regular	Black	Beige
Small	1401-UC-BKR	1411-UC-BKR
Medium	1402-UC-BKR	1412-UC-BKR
Large	1403-UC-BKR	1413-UC-BKR
X Large	1404-UC-BKR	1414-UC-BKR
XX Large	1405-UC-BKR	1415-UC-BKR
Size Tall		
Small	1401-UC-BKT	1411-UC-BKT
Medium	1402-UC-BKT	1412-UC-BKT
Large	1403-UC-BKT	1413-UC-BKT
X Large	1404-UC-BKT	1414-UC-BKT
XX Large	1405-UC-BKT	1415-UC-BKT

PRODUCT CODE		
Size Regular	Black	Beige
Small	1401-FC-BKR	1411-FC-BKR
Medium	1402-FC-BKR	1412-FC-BKR
Large	1403-FC-BKR	1413-FC-BKR
X Large	1404-FC-BKR	1414-FC-BKR
XX Large	1405-FC-BKR	1415-FC-BKR
Size Tall		
Small	1401-FC-BKT	1411-FC-BKT
Medium	1402-FC-BKT	1412-FC-BKT
Large	1403-FC-BKT	1413-FC-BKT
X Large	1404-FC-BKT	1414-FC-BKT
XX Large	1405-FC-BKT	1415-FC-BKT

See page 84 for Sizing Chart

Liner Options



Basic Liner
Non-compressive

- Non-compressive liner for the protection of the skin
- Easy to apply knee-high liners
- drirelease moisture wicking and odour control
- Flat seams for added comfort
- Silver ionized treatment for bacteriostatic properties
- Available in Heathered Silver colour only

Size: One size fits all	
Pair	10 Pair Pack
1700 - NK1	1700 - NKC
Included With Compreflex Standard Calf & Foot Compreflex Standard Calf	



Transition Liner
10–15 mmHg

- Mild compression in foot and ankle only
- drirelease moisture wicking and odour control
- Extra thick sole from toe to heel for added comfort
- Silver ionized treatment for bacteriostatic properties
- Available in Black and Beige


PRODUCT CODE (Included with Compreflex Transition)		
Size	Pair (Black)	Pair (Beige)
Small	1601-CK	1611-CK
Medium	1602-CK	1612-CK
Large	1603-CK	1613-CK
X Large	1604-CK	1614-CK
XX Large	1605-CK	1615-CK



Complete Liner
10–15 mmHg

- Mild compression in foot and ankle with graduated compression in calf
- Comfortable knee-high liners
- Drirelease moisture wicking and odour control
- Extra thick sole from toe to heel for added comfort
- Silver ionized treatment for bacteriostatic properties
- Available in Black and Beige

PRODUCT CODE (Included with Compreflex Complete)		
Size	Pair (Black)	Pair (Beige)
Small	1601-CK-F1	1611-CK-F1
Medium	1602-CK-F1	1612-CK-F1
Large	1603-CK-F1	1613-CK-F1
X Large	1604-CK-F1	1614-CK-F1
XX Large	1605-CK-F1	1615-CK-F1



Foot Size:
I: Measure circumference of ankle bend and heel
J: Measure circumference across metatarsal heads
K: Measure **medial** length from heel to 1st metatarsal head

LINER

TRANSITION AND COMPLETE LINERS SIZING CHART				
Size	Ankle Circ. (A)	Calf Circ. (C)	Women's Shoe	Men's Shoe
Small	19–26 cm	29–39 cm	<7½	<6½
Medium	21–30 cm	34–44 cm	7.5–9.5	6.5–8.5
Large	26–36 cm	40–50 cm	10–12	9–11
X Large	31–41 cm	46–56 cm	12.5–15	11.5–14
XX Large	36–46 cm	54–64 cm	12.5–15	11.5–14

Compreboot Options



Compreboot Lite **FOOT**
Variable Compression

- Ultra-thin profile that easily fits into most footwear
- Highly adjustable design may be trimmed to fit
- Provides mild, therapeutic compression

PRODUCT CODE (Available in Black Only)	
Size	Black
Small / Medium / Large	
Regular	1803-BLR
Long	1803-BLL
X Large / XX Large	
Regular	1805-BLR
Long	1805-BLL



Compreboot Standard **FOOT**
Variable Compression

- Four-way stretch for firm compression around ankle
- Ideal for mild to moderate chronic edema
- Low stretch material conforms to foot shape

PRODUCT CODE		
Size	Black	Beige
Small		
Regular	1801-BTR	1811-BTR
Long	1801-BTL	1811-BTL
Medium / Large		
Regular	1803-BTR	1813-BTR
Long	1803-BTL	1813-BTL
X Large / XX Large		
Regular	1805-BTR	1815-BTR
Long	1805-BTL	1815-BTL



Compreboot Plus **FOOT**
Variable Compression

- Chip foam lining creates high / low pressure ridges
- Removable gripper patches provide traction
- Provides moderate, therapeutic compression

PRODUCT CODE (Available in Black Only)	
Size	Black
Small	
Regular	1801-BPR
Long	1801-BPL
Medium / Large	
Regular	1803-BPR
Long	1803-BPL
X Large / XX Large	
Regular	1805-BPR
Long	1805-BPL

COMPREBOOT LITE

FOOT SIZING CHART (I, J, K)				
	Small / Medium / Large		X Large / XX Large	
	Regular	Long	Regular	Long
I	28–36 cm	28–36 cm	≤44 cm	≤44 cm
J	22–26 cm	22–26 cm	≤33 cm	≤33 cm
K	14–18 cm	19–23 cm	20–22 cm	25–27 cm

COMPREBOOT STANDARD AND PLUS

FOOT SIZING CHART (I, J, K)						
	Small		Medium / Large		X Large / XX Large	
	Regular	Long	Regular	Long	Regular	Long
I	28–36 cm	28–36 cm	<39 cm	<39 cm	<44 cm	<44 cm
J	22–26 cm	22–26 cm	<30 cm	<30 cm	<33 cm	<33 cm
K	14–18 cm	19–23 cm	18–20 cm	23–25 cm	20–22 cm	25–27 cm



Compreflex Strap Extender

- Provides up to 10 cm of additional circumferential size range to Compreflex calf wraps
- Detaches easily
- Potential solution for patients outside of standard size ranges of Compreflex calf wraps or needing custom size

	1 Unit	10 Unit Pack
Black	1401-SE1	1401-SEB
Beige	1411-SE1	1411-SEB



Accutab Devices

The Accutab accurately delivers the level of compression prescribed by the health care provider. An easy to use device, the Accutab ensures the safe delivery of compression in doses of 20–30mmHg, 30–40mmHg, or 40–50mmHg.



Standard Thigh

20–50 mmHg



Features

- The Compreflex Standard Thigh is an inelastic compression device for the thigh. The product includes a Belt and Hip Attachment to support the product throughout daily use.
- Pair the Compreflex Standard Thigh with the Compreflex Standard Knee (sold separately) to create a complete compression system for the upper leg.

Ideal for

- Chronic Venous Disease (CVD) including, but not limited to, ulcer due to venous stasis
- Lymphedema - Active and Maintenance
- Non-specific or post surgical edema

Includes

- Compreflex Standard Thigh
- 1 Pair of non-compressive Cotton Liners per leg
- Belt and Hip Attachment
- Accutab Devices

Colours




COMPREFLEX STANDARD THIGH
SIZE CHART (D-F) & PRODUCT CODE

	Small	Medium	Large	X Large
F	48–58cm	56–66cm	64–74cm	74–84cm
E	43–53cm	51–61cm	58–68cm	68–78cm
D	38–48cm	46–56cm	53–63cm	63–73cm

BLACK					
Left	Regular	1401-TCR-L	1402-TCR-L	1403-TCR-L	1404-TCR-L
	Tall	1401-TCT-L	1402-TCT-L	1403-TCT-L	1404-TCT-L
Right	Regular	1401-TCR-R	1402-TCR-R	1403-TCR-R	1404-TCR-R
	Tall	1401-TCT-R	1402-TCT-R	1403-TCT-R	1404-TCT-R

BEIGE					
Left	Regular	1411-TCR-L	1412-TCR-L	1413-TCR-L	1414-TCR-L
	Tall	1411-TCT-L	1412-TCT-L	1413-TCT-L	1414-TCT-L
Right	Regular	1411-TCR-R	1412-TCR-R	1413-TCR-R	1414-TCR-R
	Tall	1411-TCT-R	1412-TCT-R	1413-TCT-R	1414-TCT-R





COMPREFLEX STANDARD KNEE SIZE CHART (C1-D)

	Small	Medium	Large	X Large
D	38-48 cm	46-56 cm	53-63 cm	58-73 cm
C1	29-39 cm	34-44 cm	39-49 cm	44-65 cm
Black	1401-KP	1402-KP	1403-KP	1404-KP
Beige	1411-KP	1412-KP	1413-KP	1414-KP


Standard Knee

20-50 mmHg



Features



- Provides accurate, inelastic compression levels of 20-50 mmHg
- Includes inner donning sleeve for easy application
- Product includes a Butterfly Strap for hard to compress areas



*Arrows indicate reference points of D and C1. Length measurement is not required

COMPREKNEE™

Variable compression



Standard Knee
CLOSED PATELLA




COMPREKNEE SIZE CHART (C1-D)


	Small	Medium	Large / X Large
D	38-48 cm	46-56 cm	53-63 cm
C1	29-39 cm	34-44 cm	39-49 cm
Black	1101-KP	1102-KP	1103-KP
Beige	1111-KP	1112-KP	1113-KP


Features

- Garment covers the knee, providing full coverage
- Straps provide a full-range of adjustment to create a comfortable fit
- Conforming material for a firm secure fit

Colours

Black

Beige



*Arrows indicate reference points of D and C1. Length measurement is not required

COOLFLEX™

20–50 mmHg



Wrap Details

LIGHTWEIGHT, BREATHABLE COMPRESSION WEAR

Coolflex is a lightweight compression device, optimal for individuals with fragile or broken skin. Made of a non-stretch mesh material, the Coolflex transfers heat and wicks away moisture, leaving the skin cool and dry.

The features of the closure system, along with the built-in guidelines for application, promote independence and encourage improved patient outcomes.

This product is patent pending.

MOST BREATHABLE GARMENT

The Spacer Fabric® TR is designed to be highly breathable for cool-wearing comfort

ONE-HANDED DONNING

EZ-Close Pressure System of eight tabs makes one-handed, lateral donning simple, even for patients with dexterity issues

INNOVATIVE DESIGN

Durable, metal hook closures for maximum strength and durability



Ideal for

- Chronic Venous Disease (CVD) including, but not limited to, ulcer due to venous stasis
- Well-managed Lymphedema
- Non-specific or post-surgical edema


Features

- Lightweight and breathable for cool-wearing comfort
- Easy to put on for patients with limited dexterity
- Durable metal clasps, instead of Velcro closure
- Low profile design allows boot to fit in most footwear
- Spacer Fabric TR is highly breathable for cool-wearing comfort



ACCURATE COMPRESSION IS EASY WITH THE COOLFLEX

Each Coolflex calf garment includes the unique, EZ-Close™ pressure system. The EZ-Close pressure system is positioned laterally, allowing the patient to apply the product without bending forward. The bands of the closure system are made to an optimal length, making them easy to grip and close using minimal hand strength and little dexterity.



The integrated EZ-Close Pressure System has visual indicators for the accurate delivery of 20–30 mmHg, 30–40 mmHg or 40–50 mmHg on each tab



Standard Calf

20–50 mmHg



Ideal for both venous and lymphatic limb disorders, the garments are available in standard and max sizing and can accommodate limbs that feature disproportionately small ankles with a larger calf.

Features

- Lightweight and breathable for cool-wearing comfort
- Easy to put on for patients with limited dexterity
- Durable metal clasps, instead of velcro closure

Includes

- Coolflex Standard Calf
- 1 Pair of non-compressive basic Liners per leg

Colour



Black



Standard Foot


Variable compression



Foot compression garments are pictured with the Basic Liner. Basic liner is included with the wrap but not with the footpiece

Features

- Low profile design allows boot to fit in most footwear
- Spacer Fabric TR is highly breathable for cool-wearing comfort
- Durable metal clasps, instead of velcro closure




Measuring Instructions

Foot Size:

- I:** Measure circumference of ankle bend and heel
- J:** Measure circumference across metatarsal heads
- K:** Measure medial length from heel to 1st metatarsal head

COOLFLEX STANDARD CALF SIZE CHART (A-C1)						
(Available in Black Only)						
	Small	Small Max	Medium / Large	Medium / Large Max	X Large / XX Large	
C1	28–39cm	35–50cm	34–50cm	45–64cm	45–64cm	
C	28–39cm	35–50cm	34–50cm	45–64cm	45–64cm	
B	23–33cm	27–38cm	29–43cm	35–49cm	37–54cm	
A	16–26cm	16–26cm	23–33cm	23–33cm	30–42cm	
Left	Regular	2101-NF-BKR-L	2101-NF-BMR-L	2103-NF-BKR-L	2103-NF-BMR-L	2105-NF-BKR-L
	Tall	2101-NF-BKT-L	2101-NF-BMT-L	2103-NF-BKT-L	2103-NF-BMT-L	2105-NF-BKT-L
Right	Regular	2101-NF-BKR-R	2101-NF-BMR-R	2103-NF-BKR-R	2103-NF-BMR-R	2105-NF-BKR-R
	Tall	2101-NF-BKT-R	2101-NF-BMT-R	2103-NF-BKT-R	2103-NF-BMT-R	2105-NF-BKT-R



COOLFLEX STANDARD FOOT SIZE CHART (I, J, K)						
(Available in Black Only)						
	Small		Medium / Large		X Large / XX Large	
	Regular	Long	Regular	Long	Regular	Long
I	28–39cm	N/A	30–43cm	N/A	34–49cm	N/A
J	22–26cm	N/A	25–30cm	N/A	28–33cm	N/A
K	14–18cm	19–21cm	18–20cm	20–22cm	20–22cm	23–25cm
Right	2101-BTR-R	2101-BTL-R	2103-BTR-R	2103-BTL-R	2105-BTR-R	2105-BTL-R
Left	2101-BTR-L	2101-BTL-L	2103-BTR-L	2103-BTL-L	2105-BTR-L	2105-BTL-L

COMPRESHORT®

10–15mmHg



IDEAL FOR PATIENTS WITH ABDOMINAL OR TRUNCAL EDEMA WHO NEED MILD, SUPPORTIVE COMPRESSION

Product available in short and capri.



Features

The Compreshort offers mild, supportive compression for the treatment of abdominal or truncal and upper thigh edema. Available in short and capri lengths, patients with lymphedema and lipedema benefit from the therapeutic compression and support the garment provides.

Ideal for

- Lymphedema
- Lipedema


Colour



Black

COMPRESHORT SIZING CHART (Waist, Hip, Midthigh)				PRODUCT CODE	
	Waist	Hip	Midthigh	Product SKU (Short)	Product SKU (Capri)
Small	69–76cm	84–91cm	41–46cm	1301-CS	1301-CP
Small Max	76–86cm	91–102cm	41–46cm	1301-CM	1301-CPM
Medium	76–86cm	91–102cm	46–53cm	1302-CS	1302-CP
Medium Max	86–97cm	102–112cm	46–53cm	1302-CM	1302-CPM
Large	86–97cm	102–112cm	53–61cm	1303-CS	1303-CP
Large Max	97–112cm	112–127cm	53–61cm	1303-CM	1303-CPM
X Large	97–112cm	112–127cm	61–71cm	1304-CS	1304-CP
X Large Max	112–132cm	127–147cm	61–71cm	1304-CM	1304-CPM
2X Large	112–132cm	127–147cm	71–81cm	1305-CS	1305-CP
2X Large Max	132–152cm	147–168cm	71–81cm	1305-CM	1305-CPM
3X Large	132–152cm	147–168cm	81–91cm	1306-CS	1306-CP
3X Large Max	152–172cm	168–188cm	81–91cm	1306-CM	1306-CPM
4X Large	152–172cm	168–188cm	91–101cm	1307-CS	1307-CP
4X Large Max	172–193cm	188–208cm	91–101cm	1307-CM	1307-CPM
5X Large	172–193cm	188–208cm	101–112cm	1308-CS	1308-CP

Note: Available in Black Only





Help Yourself

Why choose **Accessories**?

- Aid in putting on and removing stockings
- Stocking cleaning to maximize usage
- Extra liners for added comfort

ACCESSORIES	
Product	Page
Gloves	100
Cotton Underliner	100
Sim-Slide	101
Magnide On/Off	101
Washing Solution	102
Sigvaris Doff N' Donner	102
Sigvaris Doff N' Donner Cone	102
Melany	102
It Stays	102

GLOVES



Rubber Gloves*

Ridged - The practical and protective way to apply and remove your compression stockings.



Latex-Free Gloves

Gripping donning gloves recommended for use in a latex-free environment.



Textile Gloves*

These high-quality gloves have a latex surface to provide a good grip when putting on and taking off compression stockings.

PACKAGE: 12 PER PACK	
SIZE	PRODUCT CODE
Small	591R400S
Medium	591R400M
Large	591R400L
X-Large	591R400X
ASST'D 4X S/M/L	591R400A

PACKAGE: 12 PER PACK	
SIZE	PRODUCT CODE
Small	592R400S
Medium	592R400M
Large	592R400L
X-Large	592R400X

PACKAGE: 12 PER PACK	
SIZE	PRODUCT CODE
X-Small	593R400XS
Small	593R400S
Medium	593R400M
Large	593R400L
X-Large	593R400X
ASST'D 4x XS/S/M	593R400A

RUBBER GLOVE AND LATEX-FREE GLOVE SIZING CHART			
Measure Circumference of the Hand at the Palm			
16.5–18 cm (6.5"–7")	19–20.5 cm (7.5"–8")	21.5–23 cm (8.5"–9")	24–25.5 cm (9.5"–10")
Small	Medium	Large	X-Large

TEXTILE GLOVE SIZING CHART					
Measure Length and Width of the Hand					
Length	20.5–21 cm	21.5–22 cm	22.5–23 cm	23.5–24 cm	24.5–25 cm
Width	7–7.5 cm	7.5–8 cm	8–8.5 cm	9–9.5 cm	10–10.5 cm
	X-Small	Small	Medium	Large	X-Large

LINERS



Cotton Underliner

10–15 mmHg
Can be used as an underlayer to protect fragile skin, secure wound dressings, add 10–15 mmHg of compression, or create an easy glide surface to aid in putting on and taking off compression stockings.

For sizing, please refer to the Natural Rubber Calf sizing on page 49.

PACKAGE: 2 PER PACK		
SIZE	PRODUCT CODE	NATURAL RUBBER SIZING
Small / Short	588USS	S1/S3
Small / Long	588USL	S2/S4
Medium / Short	588UMS	M1/M3
Medium / Long	588UML	M2/M4
Large-XL / Short	588ULS	L1/L3 - X1/X3
Large-XL / Long	588ULL	L2/L4 - X2/X4

DONNING AIDS



Scan here to see how to use the Sim-Slide®
bit.ly/SIGssyt



Scan here to see how to use the Magnide On/Off®
bit.ly/SIGmofyt

Sim-Slide (For open toe Styles)

SIZE	PRODUCT CODE
Small	587SLIDESM
Medium	587SLIDEMD
Large	587SLIDELG
X-Large	587SLIDEXL

Sizing Chart
For Sim-Slide

Women's Shoe Size	≤4–6	6.5–8	8.5–9.5	≥10
Men's Shoe Size	≤3–4.5	5–6.5	7–8	≥8.5
Size	S	M	L	XL

Magnide On/Off (For closed toe Styles)

SIZE	PRODUCT CODE
Medium	587MAGMD
Large	587MAGLG
X-Large	587MAGXL

SIZING CHART
For Magnide On/Off

Women's Shoe Size	≤8.5	9–12	≥ 12.5
Men's Shoe Size	≤7.5	8–11	≥ 11.5
Size	M	L	XL

* XL suitable for extra-wide feet

* Caution: Rubber Gloves and Textile Gloves contain natural rubber latex that may cause allergic reaction

DONNING AIDS



SIGVARIS Doff N' Donner

A soft device that allows you to quickly guide a compression sock or stocking up and down your leg or arm, making it simple to put on and remove all styles of graduated compression. It is perfect for both self-donning and caregiver use. Latex free.

SIZE	PRODUCT CODE
One Size	583D100



Scan here to see how to use
Sigvaris DOFF N' DONNER
bit.ly/sigdnd



SIGVARIS Doff N' Donner Cone

This loading device makes it easy to roll stockings onto the SIGVARIS Doff N' Donner. The bottom of the base firmly attaches to most clean, flat surfaces via suction.

SIZE	PRODUCT CODE
One Size	583D200



Scan here to see how to use
the MELANY
goto.sigvaris.com/bbiMdz



MELANY

MELANY from Sigvaris makes both donning and doffing compression garments simple. Its innovative one-piece design with adjustable handles was made especially for users with limited range of motion and finger dexterity.

SIZE	CALF CIRCUMFERENCE	PRODUCT CODE
X-Small	25–30 cm	579SXS
Small	30–38 cm	579SSM
Medium	38–46 cm	579SMD
Large	46–56 cm	579SLG

WASHING SOLUTION



Washing Solution

Gentle washing solution specifically formulated to help maintain the quality of SIGVARIS GROUP compression socks and stockings. Also great for washing delicate garments!

PACKAGE OPTIONS	
CONTAINS	PRODUCT CODE
Six 16 oz. Bottles (473ml)	586W900
One 16 oz. Bottle (473ml)	586W950
Twelve 2 oz. Bottles (60ml)	586W975

IT STAYS



Roll-on body adhesive

A gentle roll-on body adhesive that washes off with water and may be used to hold compression garments in place against the skin. IT STAYS! is pliable and moves with the body.

SIZE	PRODUCT CODE
One 2 oz Bottle (60ml)	581I300
Twelve 2 oz. Bottles (60ml)	581I350

SYMPTOMS / INDICATIONS	COMPRESSION LEVEL AT THE ANKLE
<ul style="list-style-type: none">Tired, achy legsOccupational or evening edemaLeg discomfort from long hours of standing or sittingPrevent varicosities & swelling during pregnancyPreventative for legs predisposed to riskReduce swelling during travelPrevent DVT during travel	15–20 mmHg
<ul style="list-style-type: none">Heavy, fatigued, aching legsMild edema in lower extremitiesMild varicosities with minimal edemaMild varicosities during pregnancyLymphedema upper extremitiesPost-procedure of small veinsPrevention or management of DVT <p><i>Note: A prescription is recommended for stockings with the indications of 20 mmHg or higher</i></p>	20–30 mmHg
<ul style="list-style-type: none">Moderate varicosities with mild to moderate edema (also during pregnancy)Post-procedure of larger veins to maintain treatment successCEAP C3 (edema), C4 (skin changes without ulceration), C5 (skin changes with healed ulceration), C6 (skin changes with active ulcerations)Management and treatment of DVT or Post-Thrombotic SyndromeSuperficial phlebitisLymphedema after decongestant therapy to maintain reductionOrthostatic / postural hypotensionPost-fracture, post-traumatic edema	30–40 mmHg
<ul style="list-style-type: none">Lymphedema after decongestant therapy to maintain reductionSevere Post-Thrombotic SyndromeSevere varicosities and/or edema	40 mmHg or higher
<ul style="list-style-type: none">Severe arterial insufficiency, intermittent claudication, ischemiaUncontrolled congestive heart failure	Contraindications
<ul style="list-style-type: none">Skin sensitivities or allergiesAdvanced neuropathy with sensory loss in leg or footDiabetes with arterial compromiseConfinement to bed or non-ambulatory use unless otherwise prescribed by the physicianEczema, hypodermatitis, dermatosclerosis, dermatitis	Cautions
<p>The therapeutic effect of graduated compression stockings is based on the compression level. With the proper CEAP classification, the correct compression level at the ankle can be chosen. For more severe symptoms, consult with a vascular specialist.</p>	
<p><i>No liability accepted for non-observance of contraindications and cautions.</i> Fronek HS, Bergan JJ, et al. The Fundamentals of Phlebology: Venous Disease for Clinicians. 2004. pg 151. Blaettner, et al. Leg Symptoms of Healthy People and Their Treatment with Compression Hosiery. Phlebology 2008. Compression therapy of the extremities. Partsch, Rabe, Stemmer. ISBN: 2.85480.770.7.</p>	

Sensitive Skin

Thermoregulating

Breathable

Moisture Wicking

Soft

Durable

Cotton

Wool

Microfiber

Natural Rubber

Largest Size Range

Strong and Durable

Sheer

Medium Sheer

Opaque

Stylish

Stylish

Elegant

Colours

Patterns

Sensinnov Grip-top Band

Anatomically Contoured

Cushioned Sole

Anti-bacterial

Hypoallergenic Natural fibers

Moisture Management

Term	Definition	Sigvaris Usage
CHRONIC VENOUS DISEASE (CVD)	A disease state comprised of a wide range of disorders including tired, achy legs, varicose veins, Chronic Venous Disease (CVD), skin issues, and venous leg ulcers (open and healed).	SIGVARIS GROUP manufactures MEDICAL and CONSUMER HEALTH CARE products to help reduce the risk of, improve symptoms of, and provide treatment for CVD.
DOUBLE-COVERED YARN	Spandex fibers are covered twice with yarn to prevent bare spandex from touching the skin for improved wearing comfort.	Used in nearly all SIGVARIS GROUP graduated compression socks and hosiery to create products that are easy to put on and take off, while also providing greater comfort.
DEEP VEIN THROMBOSIS (DVT)	A blood clot in a deep vein, most commonly in the lower leg or thigh.	SIGVARIS GROUP cares about DVT prevention. Clinical evidence supports the use of graduated compression to help prevent DVTs and the more serious complication of PE (Pulmonary Embolism), which can be life-threatening.
GRADUATED COMPRESSION VS. ANTI-EMBOLISM	Graduated compression is for ambulatory patients, whereas anti-embolism stockings are for non-ambulatory patients.	SIGVARIS GROUP offers true graduated compression, meaning products are tightest at the ankle and decrease in pressure going up the leg to help improve circulation.
HYPOALLERGENIC	Designed to minimize the chance of an allergic reaction.	SIGVARIS GROUP offers Latex free products and cotton fabrics for those with sensitive skin and allergies. Also, many SIGVARIS GROUP accessories are designed for consumers with sensitive skin.
ISO CERTIFIED	SIGVARIS GROUP is ISO (International Organization for Standardization) certified, indicating we have the required quality systems in place.	SIGVARIS GROUP follows ISO Standard 13485:2016 for our products manufactured in North and South America.
MICROFIBER	An extremely fine synthetic fiber knitted into textiles with the texture and drape of natural-fiber fabric. This fabric has excellent breathability.	Used to create products with a soft and comfortable feel. Find this fiber in 840 Soft Opaque, 830, 143 & 183 Microfiber Patterns, and 820 Microfiber.
mmHg (MILLIMETERS OF MERCURY)	Defined as the extra pressure generated by a column of mercury one millimeter high; the same unit of measurement used for blood pressure.	SIGVARIS GROUP products are offered in several ranges of compression levels based on this measurement.
MERINO WOOL	Merino is the softest and finest wool available in the world and is inherently itch-free.	SIGVARIS GROUP uses this premium wool yarn to create soft products that are thermoregulating for use in all seasons. Used in 422 Thermoregulating Wool, 421 MERINO OUTDOOR SOCKS, 240 Merino Wool, and 152 & 192 ALL-SEASON MERINO WOOL.
SEA ISLAND COTTON	Sea Island Cotton is the rarest cotton in the world and represents the finest in cotton softness.	SIGVARIS GROUP uses the highest-quality yarns in the industry to create products that are best in class. Look for this fine quality cotton in our 220 Sea Island and 151 & 191 SEA ISLAND COTTON SOCKS.
SENSINNOV	A patented SIGVARIS silicone top band that gently grips the skin and keeps the band in place all day.	This top band is featured in 780 Sheer, 840 Soft Opaque, 710 Patterns, and 820 Microfiber to help keep thigh-high stockings in place all day long.
SUPIMA COTTON	Supima accounts for only about three percent of annual cotton production in the United States. Its fineness and longer staple length make it a premium cotton fiber.	SIGVARIS GROUP uses this premium cotton yarn in the 230 Cotton, 220 Sea Island and 151 & 191 SEA ISLAND COTTON SOCKS, 146 & 186 CASUAL COTTON, and 360, 142 & 182 Cushioned Cotton products, making them comfortable and ideal for those with sensitive skin.

Please follow these instructions carefully to determine the correct size for different product styles.

Measurements should always be taken in the following order, and as early in the morning as possible when swelling is minimal.

CALF STYLE
Open/Closed Toe

1	Measure the ankle circumference at the narrowest part of the ankle, above the ankle bone. Use the appropriate series sizing chart to determine if ankle is SMALL, MEDIUM, LARGE or X-LARGE. Note: This should be your first measurement.
2	Measure the calf circumference at the fullest part of the calf.
3	Determine the calf length by measuring from the floor up to the bend of the knee. This measurement is taken on the outside of the leg.



1 Ankle circumference directly above ankle bone



2 Calf circumference at fullest part of the calf

PANTYHOSE AND THIGH-HIGH STYLE
Open/Closed Toe

1	Measure the ankle circumference at the narrowest part of the ankle, above the ankle bone. Use the appropriate series sizing chart to determine if ankle is SMALL, MEDIUM or LARGE. Note: This should be your first measurement.
2	Measure the calf circumference at the fullest part of the calf.
3	Skip Step 3 (but refer to General Guidelines #6 below).
4	Measure the thigh at its widest circumference.
5	Determine the leg length by measuring from the floor up to the gluteal fold.
6	Determine the hip circumference by measuring the widest part of the hips.



3 Calf length from the floor to the right-angle bend of the knee



4 Widest circumference of the thigh

GENERAL GUIDELINES

1	Measurements should always be taken as early in the morning as possible, when swelling is minimal.
2	Measure against bare skin, not over clothing (Steps 1-5).
3	The measuring system is based on having the foot flat against the floor and the calf or leg at a right angle to the foot (see images to right for examples).
4	Measurements should be taken by pulling the tape around the leg so that it is snug but not constricting.
5	Length measurements should be taken from the ground up.
6	Take all measurements, even if the patient only requires calf-length products, as the patient might decide at a later date to purchase additional styles.



5 Leg length from the floor up to the gluteal fold



6 Circumference at the hip

HOW TO ORDER MEDICAL AND SPECIALTY PRODUCTS

Series Number	Compression Level	Product Style	Size	Gender	Colour
23 = Cotton (See individual pages for series number)	1 = 15-20 mmHg 2 = 20-30 mmHg 3 = 30-40 mmHg 4 = 40-50 mmHg* 5 = 50-60 mmHg*	C = Knee-high or Calf (See individual pages for available styles)	ML = Medium Long (See individual pages for sizing)	M = Men W = Women O = Open toe	99 = Black (See individual pages for available colours)
23	2	C	ML	M	99

* 40-50 mmHg and 50-60 mmHg available in Natural Rubber only

To order Plus Calf with Grip-top or Plus Size Pantyhose (860P for women only):
add **"/P"** after the colour in the ordering code

To order Calf Style with Grip-top:
add **"/S"** after the colour in the ordering code (see product pages for availability)

To order Thigh with Waist Attachment:
add **"/L"** to indicate the Left leg or **"/R"** to indicate the Right leg after the colour in the ordering code (860W and 500W only)

To order Motion High Tech and Motion Thermoregulating Wool:
The letter for Gender is not used in the ordering code.
Example: 412CSM59 = Motion High Tech, 20-30 mmHg, Closed toe Calf, Size Small Medium, Colour Red

HOW TO ORDER CONSUMER HEALTH PRODUCTS

Series Number	Product Style	Size	Colour
120 = SHEER FASHION (See individual pages for series number)	N = Thigh-high (See individual pages for available styles)	A = Small (See individual pages for sizing)	33 = Honey (See individual pages for available colours)
120	N	A	33

To order ATHLETIC RECOVERY SOCKS, PERFORMANCE SLEEVES, MERINO OUTDOOR SOCKS, and EVERSOFT DIABETIC SOCKS:

Size is only one letter and the letter for Gender is not used in the ordering code.

Example: 421CL31 = MERINO OUTDOOR SOCK, 15-20 mmHg, Closed toe Calf, Size Large, Colour Olive

Sigvaris, SIGVARIS GROUP, PFS², Sensinno, Compreflex, Compreboot, and Accutab are registered trademarks of SIGVARIS AG, CH-9014 St.Gallen/Switzerland, in many countries worldwide. © Copyright 2021 SIGVARIS, Inc.

Coolflex, Spacer Fabric TR, and CoverOne are trademarks of SIGVARIS AG, CH-9014, St. Gallen/Switzerland, in many countries worldwide.

Lycra is a registered trademark of INVISTA. Supima is a registered trademark of the Supima Association of America Corporation. NILIT is a registered trademark of Nilit. All rights reserved. drirelease is copyright ©2021 drirelease. drirelease is an Optimizer Brands company. All rights reserved. Sim-Slide and Magnide are registered trademarks of Arion Holding BV.

SIGVARIS GROUP

Scan to see
digital catalog



To Order :

1.800.363.4999

Fax: 1.800.263.8736

canadaorders@sigvaris.com

Online Dealer Portal: sigvaris-online.com



sigvarisgroup.canada



[sigvarisgroup_canada](https://www.instagram.com/sigvarisgroup_canada)



[SIGVARIS GROUP Canada](https://www.linkedin.com/company/SIGVARIS%20GROUP%20Canada)



[SIGVARIS GROUP Canada](https://www.youtube.com/SIGVARIS%20GROUP%20Canada)