

SIGVARIS





SIGVARIS GROUP

2019 / 2020
COMPRESSION SOLUTIONS



Life changes. Your lifestyle doesn't have to.

At Sigvaris, we care about people. For 155 years, our compression products have offered effective and reliable solutions for venous conditions – anywhere, anytime, and for anyone.

As a Sigvaris Certified Fitter, when you provide **the Sigvaris solution** to your clients, you're helping them move with life, without compromise. Because you care about people too. And that's why we're proud to work with Fitters like you.

Our shared passion for helping people is what got us here, and we're excited to take our next step towards wellness together.

Let's take a look at what we've got in store.

Live your life.

You know your clients. They trust you to help them feel their best.

Sigvaris understands that people want to be their best selves – everywhere, every day. So we're making it even easier for you to find the right solutions for everyone.

You told us what your clients need and we are proud to continue delivering the **FREEDOM, SECURITY, and SELF-CONFIDENCE** that keeps them coming back.

Our products aren't changing. We are still the forward-thinking leader you rely on and trust. You know that Sigvaris Medical products are made in our own factories, from premium double-covered spandex yarns and fibres. Our sizing system ensures exceptional comfort, and we guarantee the compression of our socks for six months – even after they are machine washed and dried. But what your clients don't know about us might surprise them.

The Sigvaris Solution:

- Medical effectiveness and reliable relief.
- Best-in-class quality, craftsmanship, and durability.
- Premium yarns like Merino wool and Sea Island cotton.
- Freedom of movement, comfort, and fit.

Previous	Series Number	New	Segment
COTTON	230	Cotton	Essential
SELECT COMFORT	860	Opaque	
SEA ISLAND COTTON	220	Sea Island Cotton	Style
ALL-SEASON MERINO WOOL	240	Merino Wool	
ALLURE	710	Patterns	
MIDSHEER	750	Medium Sheer	
EVERSHEER	780	Sheer	
MIDTOWN MICROFIBER	820	Microfiber	
MICROFIBER SHADES	830	Microfiber Patterns	
SOFT OPAQUE	840	Soft Opaque	
CUSHIONED COTTON	360	Cushioned Cotton	Motion
TRAVERSE SOCKS	412	High Tech	
MERINO WOOL	422	Thermoregulating Wool	
SWISS MOUNTAIN SOCK	462	Swiss Mountain Sock	
DAILY COMFORT	850	Comfort	

Improved shopability



NEW PACKAGE

FRONT



BACK





Be yourself.

ESSENTIAL

Support Yourself

Sigvaris Essential solutions always meet your clients' needs. These straightforward stockings are an everyday go-to for people who seek the comfort, support and reliability of traditional therapeutic socks and stockings.

Why choose **Essential**?

- Reduce swelling and discomfort.
- Long-distance travel.
- Everyday comfort and reliability.

STYLE

Express Yourself

Choice is empowering, and fashion can be as unique as you are. With Sigvaris Style, personality will shine as your clients select from distinctive hosiery and socks.

Why choose **Style**?

- Fashionable, modern patterns and designs.
- Special occasion and formal wear.
- Our most luxurious feel.

MOTION

Push Yourself

When daily comfort for hard-working legs is required, put your clients in motion. Sigvaris Motion socks offer all day strong support for people with active lives.

Why choose **Motion**?

- Cushioned sole for added comfort.
- Durable thermoregulating wool.
- Moisture-wicking yarns.

MEDICAL COMPRESSION

From 15–20mmHg to 50–60mmHg

Essential

06

230 Cotton

20–30mmHg / 30–40mmHg

Closed toe: calf, calf w/grip-top, thigh-high w/grip-top

Open toe: calf, thigh-high w/grip-top

08

860 Opaque

20–30mmHg / 30–40mmHg

Closed toe: calf, petite calf, calf plus, calf w/grip-top, thigh w/grip-top, petite thigh-high w/grip-top, extra-large thigh w/grip-top, pantyhose, petite pantyhose, plus-sized pantyhose

Open toe: calf, calf plus, calf w/grip-top, thigh w/grip-top, extra-large thigh w/grip-top, thigh-high w/waist attachment

Style

12

780 Sheer

15–20mmHg / 20–30mmHg / 30–40mmHg

Closed toe: calf, thigh-high w/grip-top, pantyhose, maternity pantyhose

Open toe: calf, thigh-high w/grip-top, pantyhose

14

750 Medium Sheer

20–30mmHg

Closed toe: calf, thigh-high w/grip-top, pantyhose

16

840 Soft Opaque

15–20mmHg / 20–30mmHg / 30–40mmHg

Closed toe: calf, calf w/grip-top, thigh-high w/grip-top, pantyhose, maternity pantyhose

Open toe: calf, thigh-high w/grip-top, pantyhose

18

710 Patterns

15–20mmHg / 20–30mmHg

Closed toe: thigh-high w/grip-top, pantyhose

20

220 Sea Island Cotton

20–30mmHg

Closed toe: calf

22

240 Merino Wool

20–30mmHg

Closed toe: calf

24

820 Microfiber

15–20mmHg / 20–30mmHg / 30–40mmHg

Closed toe: calf

26

830 Microfiber Patterns

20–30mmHg

Closed toe: calf

Motion

30

360 Cushioned Cotton

20–30mmHg

Closed toe: calf

32

412 High Tech

20–30mmHg

Closed toe: calf

34

422 Thermoregulating Wool

20–30mmHg

Closed toe: calf

36

462 Swiss Mountain Sock

20–30mmHg

Closed toe: calf

38

850 Comfort

20–30mmHg

Closed toe: calf

Specialty

42

500 Natural Rubber

30–40mmHg / 40–50mmHg / 50–60mmHg

Open toe: calf, thigh-high, thigh-high w/grip-top, thighw/ waist attachment, pantyhose

44

910 Advance Armsleeve

15–20mmHg / 20–30mmHg / 30–40mmHg

Armsleeve: grip-top w/out gauntlet, grip-top w/gauntlet

76

Inelastic Wraps

84

Accessories

87

Symptoms & Indications

88

Glossary of Terms

89

Glossary of Icons

90

How to Measure

91

How to Order

92

Why SIGVARIS GROUP

CONSUMER HEALTH CARE COMPRESSION PRODUCTS

From 8–15mmHg to 20–30mmHg

48

146/186 CASUAL COTTON

15–20mmHg **Closed toe:** calf

50

142/182 CUSHIONED COTTON

15–20mmHg **Closed toe:** calf

52

151/191 SEA ISLAND COTTON

15–20mmHg **Closed toe:** calf

54

152/192 ALL-SEASON MERINO WOOL

15–20mmHg **Closed toe:** calf

56

143/183 MICROFIBER SHADES

15–20mmHg **Closed toe:** calf

58

120 SHEER FASHION

15–20mmHg **Closed toe:** calf, thigh-high w/grip-top, pantyhose, maternity pantyhose

Open toe: calf

60

170 SOFT SILHOUETTE LEGGINGS

15–20mmHg **Open foot:** leggings

62

189 BUSINESS CASUAL

15–20mmHg **Closed toe:** calf

64

401 ATHLETIC RECOVERY SOCKS

15–20mmHg **Closed toe:** calf

66

421 MERINO OUTDOOR SOCKS

15–20mmHg **Closed toe:** calf

68

160 EVERSOFIT DIABETIC SOCKS

8–15mmHg **Closed toe:** calf

70

412 PERFORMANCE SLEEVES

20–30mmHg **Open foot:** sleeve

72

TRAVENO TRAVEL SOCKS

15–20mmHg **Closed toe:** calf

sigvaris

Essential Products			
Previous	Series Number	New	Page
COTTON	230	Cotton	6
SELECT COMFORT	860	Opaque	8

Essential

ESSENTIAL

Support Yourself

Sigvaris Essential solutions always meet your clients' needs. These straightforward stockings are an everyday go-to for people who seek the comfort, support and reliability of traditional therapeutic socks and stockings.

Why choose **Essential**?

- Reduce or completely eliminate edema.
- Long-distance travel.
- Everyday comfort and reliability.

230

Essential COTTON

for women & men

Essential Cotton combines all-day cotton comfort and a classic ribbed design that is ideal for sensitive skin.



**Sensitive
Skin**



Cotton



Breathable

Compression Level

20-30mmHg

30-40mmHg

Features and Benefits

- Supima® cotton for highest level of comfort and breathability
- Ideal for sensitive skin and allergies
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off
- Soft to the skin with a silicone dot band to prevent slipping

Fibre Content

- 50% Nylon, 25% Spandex, 25% Supima Cotton (content can vary slightly by style)
- Latex free

Best Suited for

- Containment of various types of edema, including venous edema, phlebolymphedema, and early stages of lymphedema
- All-day cotton comfort in all types of climates
- Ideal for dress or casual shoes
- Ideal for long-distance travel



Silicone
dotted grip-top
band in calf
and thigh







Sizing Chart

	Small		Medium		Large		X Large	
Ankle Circumference	18–21.5cm (7"–8.5")		21.5–25.5cm (8.5"–10")		25.5–29cm (10"–11.5")		29–36cm (11.5"–14")	
Calf Circumference	28–39cm (11"–15.5")		31.5–45cm (12.5"–17.5")		35.5–51cm (14"–20")		41–61cm (16"–24")	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")
Thigh Circumference	45.5–66cm (18"–26")		48–73cm (19"–28.5")		52cm–84cm (20.5"–33")		NA	
Leg Length: Women & Open Toe	Short	Long	Short	Long	Short	Long		
	≤72cm (28.5")	>72cm (28.5")	≤72cm (28.5")	>72cm (28.5")	≤72cm (28.5")	>72cm (28.5")		
Leg Length: Men	Short	Long	Short	Long	Short	Long		
	≤76cm (30")	>76cm (30")	≤76cm (30")	>76cm (30")	≤76cm (30")	>76cm (30")		
Size	SS	SL	MS	ML	LS	LL	XS	XL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12 Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

Styles and Colours

Closed Toe = ● Open Toe = ○						
Calf (C)	○	●	●	●○	●	●
Calf w/Grip-Top (C/S) <small>Silicone grip-top band</small>				●*		●**
Thigh w/Grip-Top (N) <small>Silicone grip-top band</small>				●○		●

* 20–30mmHg only
** 30–40mmHg calf with grip-top only available in black
30–40mmHg extra large calf only available in black with grip-top

Product Code Example: 232NMLW66

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour
23	2	N	ML	W	66

Product Code Example: 232CLSM66/S

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour	Grip-top
23	2	C	LS	M	66	/S

Essential

OPAQUE

for women & men

Essential Opaque offers durability you can trust in our largest range of sizes, making it your go-to choice for everyday hosiery.



Durable



**Largest
Size Range**



Opaque

Compression Level

20-30mmHg

30-40mmHg

Features and Benefits

- Largest size range
- Durable construction combined with a smooth, opaque look
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off
- Soft to the skin with a silicone dot band to prevent slipping
- Cotton-lined gusset for pantyhose styles
- Plus size pantyhose great alternative for maternity

Fibre Content

- 64% Nylon, 36% Spandex
(content can vary slightly by style)
- Latex free

Best Suited for

- Containment of various types of edema including venous edema, phlebolympheidema, and early stages of lymphedema
- Post-vein treatment including thermal and non-thermal techniques
- Patients looking for petite and plus-sized options



Silicone
dotted grip-top
band in calf
and thigh

Sizing Chart

	Small			Medium			Large		X Large	
Ankle Circumference	18–21.5cm (7"–8.5")			21.5–25.5cm (8.5"–10")			25.5–29cm (10"–11.5")		29–36cm (11.5"–14")	
Calf Circumference	28–39cm (11"–15.5")			31.5–45cm (12.5"–17.5")			35.5–51cm (14"–20")		41–61cm (16"–24")	
Plus Size Calf Circumference	NA			39.5–53cm (15.5"–21")			45.5–61cm (18"–24")		51–71cm (20"–28")	
Calf Length	Petite*	Short	Long	Petite*	Short	Long	Short	Long	Short	Long
	≤35.5cm (14")	≤41cm (16")	>41cm (16")	≤35.5cm (14")	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")
Thigh Circumference	45.5–66cm (18"–26")			48–73cm (19"–28.5")			52–84cm (20.5"–33")		64–95cm (25"–37.5")	
Leg Length	Petite*	Short	Long	Petite*	Short	Long	Short	Long	Short	Long
	≤65cm (25.5")	≤72cm (28.5")	>72cm (28.5")	≤65cm (25.5")	≤72cm (28.5")	>72cm (28.5")	≤72cm (28.5")	>72cm (28.5")	≤72cm (28.5")	>72cm (28.5")
Hip Circumference Guidelines	83–112.5cm (32.5"–44.5")			87–127cm (34.5"–50")			98–150.5cm (38.5"–59.5")		NA	
Size	SP*	SS	SL	MP*	MS	ML	LS	LL	XS	XL
Plus Size Hip Circumference Guidelines	<142cm (56")			<157.5cm (62")			<178cm (70")		NA	
Plus Size Calf				MS/P ML/P		LS/P	LL/P	XS/P	XL/P	

Styles and Colours

Closed Toe = ● Open Toe = ○	<div>NATURAL 33</div> <div>SUNTAN 36</div> <div>CRISPA 66</div> <div>BLACK 99</div>			
Women's Calf (C)	●	●	○	●●
Men's Calf (C) (20–30mmHg only)			●	●
Calf Plus w/Grip-Top (C/P**) <small>Silicone grip-top band</small>			●○	
Calf w/Grip-Top (C/S) <small>Silicone grip-top band</small>			●○	●
Thigh w/Grip-Top (N) <small>Silicone grip-top band</small>	●	●	○	●
Size XS/XL Thigh w/Grip-Top (N) <small>Silicone grip-top band</small>			●○	
Pantyhose (P)	●	●		●
Pantyhose Plus (P/P)***			●	●
Thigh w/Waist Attachment (W) <small>Indicate leg LEFT (/L) or RIGHT (/R)</small>			○	

Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14 (Calf style only)

Women's Shoe Sizes: Petite Leg 4–8, Short Leg 5.5–9.5, Long Leg 8–12

* Petite sizing (SP & MP): available in closed toe only; without grip-top, 20–30mmHg only; Natural, Suntan, and Black only; unavailable in Plus sizes

Note: 860 closed toe calf products are sized for women and men. All other 860 closed toe products are sized for women, open toe products are unisex

** Calf plus style available only with grip-top

*** Pantyhose plus is an alternative for maternity

Product Code Example: 863PMSW33

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour
86	3	P	MS	W	33

Product Code Example: 862CLLW66/P

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour	Plus size
86	2	C	LL	W	66	/P

sigvaris

Style Products

Previous	Series Number	New	Page
EVERSHEER	780	Sheer	12
MIDSHEER	750	Medium Sheer	14
SOFT OPAQUE	840	Soft Opaque	16
ALLURE	710	Patterns	18
SEA ISLAND COTTON	220	Sea Island Cotton	20
ALL-SEASON MERINO WOOL	240	Merino Wool	22
MIDTOWN MICROFIBER	820	Microfiber	24
MICROFIBER SHADES	830	Microfiber Patterns	26

Style

STYLE

Express Yourself

Choice is empowering, and fashion can be as unique as you are. With Sigvaris Style, personality will shine as your clients select from distinctive hosiery and socks.

Why choose **Style**?

- Fashionable, modern patterns and designs.
- Special occasion and formal wear.
- Our most luxurious feel.

780

Style SHEER

for women

Style Sheer hosiery offers fashionable sheerhness you can feel comfortable in, making it ideal for dress and professional wear.



Sheer



Elegant



Breathable

Compression Level

15-20mmHg

20-30mmHg

30-40mmHg

Features and Benefits

- Fashionable for professional and dress wear
- Consistent sheerhness throughout the leg
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off
- Comfortable Sensinnov® top band prevents slipping
- Cotton-lined gusset for pantyhose styles

Fibre Content

- 15-20, 20-30mmHg: 69% Nylon, 31% Spandex (content can vary slightly by style)
- 30-40mmHg: 60% Nylon, 40% Spandex (content can vary slightly by style)
- Latex free

Best Suited for

- Management of early stages of chronic venous disease
- Pre and post vein treatment with minimal edema
- Patients who prefer fashionable options
- Pregnancy
- Ideal for business or dress



Sensinnov
grip-top band
in thigh

Sizing Chart

	Small		Medium		Large	
Ankle Circumference	18–21.5cm (7"–8.5")		21.5–25.5cm (8.5"–10")		25.5–29cm (10"–11.5")	
Calf Circumference	28–39cm (11"–15.5")		31.5–45cm (12.5"–17.5")		35.5–51cm (14"–20")	
Calf Length	Short	Long	Short	Long	Short	Long
	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")
Thigh Circumference	45.5–66cm (18"–26")		48–73cm (19"–28.5")		52–84cm (20.5"–33")	
Leg Length	Short	Long	Short	Long	Short	Long
	≤72cm (28.5")	>72cm (28.5")	≤72cm (28.5")	>72cm (28.5")	≤72cm (28.5")	>72cm (28.5")
Hip Circumference Guidelines	83–112.5cm (32.5"–44.5")		87–127cm (34.5"–50")		98–150.5cm (38.5"–59.5")	
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes: Short 5.5–9.5, Long 8–12
Maternity Hip Circumference Guidelines: Small 142cm (56"), Medium 157.5cm (62"), Large 178cm (70")

Styles and Colours

Closed Toe = ● Open Toe = ○	DARK NAVY 08	NATURAL 33	WARM SAND 34*	SUNTAN 36	TOASTED ALMOND 37*	CAFÉ 73*	MOCHA 85*	NIGHT-SHADE 94	BLACK 99
Calf (C)	●	○○	○○	○○	○○	○○	○○	●	○○
Thigh w/Grip-Top (N) <small>Sensinnov grip-top band</small>	●	○○	○○	○○	○○	○○	○○	●	○○
Pantyhose (P)	●	○○	○○	○○	○○	○○	○○	●	○○
Maternity Pantyhose (M)*		●							●

* 20–30mmHg ONLY

Product Code Example: 781NMLW36

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour
78	1	N	ML	W	36

Product Code Example: 783PSSO33

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour
78	3	P	SS	O	33

750

Style MEDIUM SHEER

for women

Thanks to the durability of Style Medium Sheer, you will have the confidence to wear sheer hosiery all day, every day.



**Medium
Sheer**



Durable



Breathable

Compression level

20-30mmHg

Features and Benefits

- Fashionable, durable, medium sheer hosiery for everyday wear
- Excellent breathability and moisture management
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off
- Soft to the skin with a silicone dot band to prevent slipping
- Cotton-lined gusset for pantyhose styles

Fibre Content

- 68% Nylon, 32% Spandex
(Content can vary slightly by style)
- Latex free

Best Suited for

- Management of early stages of chronic venous disease
- Pre and post vein treatment with minimal edema
- Daily wear for business or casual
- Active professionals needing a durable sheer stocking



Silicone
dotted grip-top
band in thigh

Sizing Chart

	Small		Medium		Large	
Ankle Circumference	18–21.5cm (7"–8.5")		21.5–25.5cm (8.5"–10")		25.5–29cm (10"–11.5")	
Calf Circumference	28–39cm (11"–15.5")		31.5–45cm (12.5"–17.5")		35.5–51cm (14"–20")	
Calf Length	Short	Long	Short	Long	Short	Long
	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")
Thigh Circumference	45.5–66cm (18"–26")		48–73cm (19"–28.5")		52–84cm (20.5"–33")	
Leg Length	Short	Long	Short	Long	Short	Long
	≤72cm (28.5")	>72cm (28.5")	≤72cm (28.5")	>72cm (28.5")	≤72cm (28.5")	>72cm (28.5")
Hip Circumference Guidelines	83–112.5cm (32.5"–44.5")		87–127cm (34.5"–50")		98–150.5cm (38.5"–59.5")	
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes: Short leg 5.5–9.5, Long leg 8–12

Styles and Colours

Closed Toe = ●				
<div><div>DARK NAVY 08</div><div>NATURAL 33</div><div>NIGHT-SHADE 94</div><div>BLACK 99</div></div>				
Calf (C)	●	●	●	●
Thigh w/Grip-Top (N) <small>Silicone grip-top band</small>	●	●	●	●
Pantyhose (P)	●	●	●	●

Product Code Example: 752NMLW33

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour
75	2	N	ML	W	33

Style **SOFT OPAQUE** for women

Style Soft Opaque hosiery blends a fashionable opaque look with long-lasting microfiber softness for your professional and everyday wear.



Opaque



Microfiber



Soft

Compression Level

15-20mmHg

20-30mmHg

30-40mmHg

Features and Benefits

- Fashionable opaque hosiery for professional and everyday wear
- Lasting microfiber softness and excellent moisture management
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off
- Comfortable Sensinnov top band on thigh style prevents slipping
- Cotton-lined gusset for pantyhose styles

Fibre Content

- 79% Nylon, 21% Spandex
(Content can vary slightly by style)
- Latex free

Best Suited for

- Management of early stages of chronic venous disease
- Pre and post vein treatment with minimal edema
- Patients who prefer fashionable options
- Pregnancy
- Long distance travel



Sensinnov
grip-top band
in thigh



Silicone
dotted grip-top
band available
in calf

Sizing Chart

	Small		Medium		Large	
Ankle Circumference	18–21.5cm (7"–8.5")		21.5–25.5cm (8.5"–10")		25.5–29cm (10"–11.5")	
Calf Circumference	28–39cm (11"–15.5")		31.5–45cm (12.5"–17.5")		35.5–51cm (14"–20")	
Calf Length	Short	Long	Short	Long	Short	Long
	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")
Thigh Circumference	45.5–66cm (18"–26")		48–73cm (19"–28.5")		52–84cm (20.5"–33")	
Hip Circumference Guidelines	45.5–66cm (32.5"–44.5")		87–127cm (34.5"–50")		98–150.5cm (38.5"–59.5")	
Leg Length	Short	Long	Short	Long	Short	Long
	≤72cm (28.5")	>72cm (28.5")	≤72cm (28.5")	>72cm (28.5")	≤72cm (28.5")	>72cm (28.5")
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes: Short leg 5.5–9.5, Long leg 8–12 Maternity hip circumference guidelines: Small: <142cm (<56"), Medium: <157.5cm (<62"), Large: <178cm (<70")

Styles and Colours

Closed Toe = ● Open Toe = ○	MIDNIGHT BLUE 09	NUDE 35	PECAN 40	TEAL 51*	MULBERRY 60*	BLUE IRIS 68*	CARAMEL 83*	ESPRESSO 89*
	●	○	●	●	●	●	●	●
Calf (C)	●	○	●	●	●	●	●	●
Calf w/Grip-Top (C/S)*								
Thigh w/Grip-Top (N) <small>Sensinnov grip-top band</small>	●	○	●	●	●	●	●	●
Pantyhose (P)	●	○	●	●	●	●	●	●
Maternity Pantyhose (M)*								

* 20–30mmHg ONLY

Closed Toe = ● Open Toe = ○	GRAPHITE 91	MIST 92*	BLACK 99
	●	●	●○
Calf (C)	●	●	●○
Calf w/Grip-Top (C/S)*			●
Thigh w/Grip-Top (N) <small>Sensinnov grip-top band</small>	●	●	●○
Pantyhose (P)	●	●	●○
Maternity Pantyhose (M)*			●

Product Code Example: 842CMSW68

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour
84	2	C	MS	W	68

Style PATTERNS

for women

Style Pattern hosiery is available in fashionable patterns, making it a staple for your professional and everyday wear.

Compression Level

15-20mmHg

20-30mmHg

Features and Benefits

- Stylish diamond-dot and stripe patterns for professional and dress wear
- Excellent breathability and moisture management
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off
- Comfortable Sensinnov top band prevents slipping
- Cotton-lined gusset for pantyhose styles

Fibre Content

- 64% Nylon, 36% Spandex
(Content can vary slightly by style)
- Latex free

Best Suited for

- Management of early stages of chronic venous disease
- Pre and post vein treatment with minimal edema
- Patients who prefer fashionable options
- Ideal for business or dress



Sheer



Elegant



Breathable



NEW!



Scalloped Sensinnov
thigh band on
Stripe patterns*





Sensinnov
thigh band on
Diamond Dot
patterns

Sizing Chart

	Small		Medium		Large	
Ankle Circumference	18–21.5cm (7"–8.5")		21.5–25.5cm (8.5"–10")		25.5–29cm (10"–11.5")	
Calf Circumference	28–39cm (11"–15.5")		31.5–45cm (12.5"–17.5")		35.5–51cm (14"–20")	
Thigh Circumference	45.5–66cm (18"–26")		48–73cm (19"–28.5")		52–84cm (20.5"–33")	
Leg Length	Short	Long	Short	Long	Short	Long
	≤72cm (28.5")	>72cm (28.5")	≤72cm (28.5")	>72cm (28.5")	≤72cm (28.5")	>72cm (28.5")
Hip Circumference Guidelines	83–112.5cm (32.5"–44.5")		87–127cm (34.5"–50")		98–150.5cm (38.5"–59.5")	
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes: Short leg 5.5–9.5, Long leg 8–12.

Styles and Colours

Closed Toe = ●					
	DARK NAVY STRIPE 13*	GRAPHITE STRIPE 17*	ONYX STRIPE 19*	GRAPHITE 91	BLACK 99
Thigh w/Grip-Top (N) Sensinnov grip-top band	●	●	●	●	●
Pantyhose (P)	●	●	●	●	●

* Coming in Fall 2019

Product Code Example: 712NLSW17

Series	mmHg Level	Style	Size	Gender/Open Toe	Colour
71	2	N	LS	W	17

220

Style

SEA ISLAND COTTON

for women & men

Style Sea Island Cotton socks are made with the world's rarest cotton, ensuring a luxurious feel that is as soft as cashmere.



Soft



Cotton



Breathable

Compression Level

20-30mmHg

Features and Benefits

- Knitted with the world's rarest cotton for a soft, luxurious feel
- Ideal for sensitive skin
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off

Fibre Content

- 43% Sea Island Cotton, 39% Nylon, 11% Spandex, 7% Supima Cotton
- Latex free

Best Suited for

- Management of symptoms related to chronic venous disease
- Patients who prefer natural fibres
- Long-distance travel



Sizing Chart

	Small		Medium		Large		X Large	
Ankle Circumference	18–21.5cm (7"–8.5")		21.5–25.5cm (8.5"–10")		25.5–29cm (10"–11.5")		29–36cm (11.5"–14")	
Calf Circumference	28–39cm (11"–15.5")		31.5–45cm (12.5"–17.5")		35.5–51cm (14"–20")		41–61cm (16"–24")	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")
Size	SS	SL	MS	ML	LS	LL	XS	XL

Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14 Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

Styles and Colours

Closed Toe = ●	<div><div>NAVY 10</div><div>BROWN 11</div><div>BLACK 99</div></div>		
	Calf (C)	●	●

Product Code Example: 222CLLM99

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour
22	2	C	LL	M	99

Style

MERINO WOOL

for women & men

Style Merino Wool socks combine extra-fine Australian Merino wool with a stylish design, making it the ultimate sock for year-round wear.



Wool



Thermoregulating



Stylish

Compression Level

20-30mmHg

Features and Benefits

- Made with the highest quality extra-fine Merino wool
- Thermoregulating for year-round comfort
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off

Fibre Content

- 53% Merino Wool, 38% Nylon, 9% Spandex
- Latex free

Best Suited for

- Management of symptoms related to chronic venous disease
- Patients who prefer fashionable and premium fabric options
- Long-distance travel



Sizing Chart

	Small		Medium		Large		X Large	
Ankle Circumference	18–21.5cm (7"–8.5")		21.5–25.5cm (8.5"–10")		25.5–29cm (10"–11.5")		29–36cm (11.5"–14")	
Calf Circumference	28–39cm (11"–15.5")		31.5–45cm (12.5"–17.5")		35.5–51cm (14"–20")		41–61cm (16"–24")	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")
Size	SS	SL	MS	ML	LS	LL	XS	XL

Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14 Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

Styles and Colours

Closed Toe = ●	<div><div>NAVY 10</div><div>BROWN 11</div><div>BLACK 99</div></div>		
	Calf (C)	●	●

Product Code Example: 242CMSW11

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour
24	2	C	MS	W	11

820

Style MICROFIBER

for men

Style Microfiber socks offer microfiber softness for your professional and everyday wear.



Microfiber



Soft



Stylish

Compression Level

15-20mmHg

20-30mmHg

30-40mmHg

Features and Benefits

- 100% of microfiber sits next to skin for a soft feel
- Excellent breathability and moisture management for added comfort
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off
- Sensinno grip-top band on thigh style

Fibre Content

- 87% Nylon , 13% Spandex
(content can vary slightly by style)
- Latex free

Best Suited for

- Management of early stages of chronic venous disease
- Patients who prefer fashionable options
- Men with larger feet
- Long-distance travel



Silicone
dotted grip-top
band in calf

Sensinno
grip-top band
in thigh

Sizing Chart

	Small		Medium		Large		X Large	
Ankle Circumference	18–21.5cm (7"–8.5")		21.5–25.5cm (8.5"–10")		25.5–29cm (10"–11.5")		29–36cm (11.5"–14")	
Calf Circumference	28–39cm (11"–15.5")		31.5–45cm (12.5"–17.5")		35.5–51cm (14"–20")		41–61cm (16"–24")	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")
Thigh Circumference	45.5–66cm (18"–26")		48–73cm (19"–28.5")		52–84cm (20.5"–33")		NA	
Leg Length	Short	Long	Short	Long	Short	Long		
	≤76cm (30")	>76cm (30")	≤76cm (30")	>76cm (30")	≤76cm (30")	>76cm (30")		
Size Shoe size	SS 7–11	SL 9.5–14	MS 7–11	ML 9.5–14	LS 7–11	LL 9.5–14	XS 7–11	XL 9.5–14
				MX* 12.5–17.5		LX* 12.5–17.5		XX* 12.5–17.5

Men's shoe sizes: Short leg 7–11, long leg 9.5–14 (For extended foot sizes: long leg 12.5–17.5)
* MX, LX, and XX extended foot sizes only available in calf style in Black & Tan Khaki in 15–20 & 20–30mmHg only. Extended foot sizes not available with grip top

Styles and Colours

Closed Toe = ●	<div><div>NAVY 10***</div><div>TAN KHAKI 32</div><div>STEEL GREY 93***</div><div>BLACK 99</div></div>			
Calf (C)	●	●	●	●
Calf w/Grip-Top (C/S)** <small>Silicone grip-top band</small>		●		●
Thigh w/Grip-Top (N) <small>Sensinnov grip-top band</small>		●		●

** Calf w/ Grip-Top not available in 15–20mmHg
*** Navy and Steel Grey colours only available in 15–20mmHg and 20–30mmHg

Product Code Example: 822NLSM99

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour
82	2	N	LS	M	99

Product Code Example: 822CMLM32/S

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour	Grip-top
82	2	C	ML	M	32	/S

Style

MICROFIBER PATTERNS

for women & men

Style Microfiber Patterns socks combine ultimate comfort with a variety of fashionable colours and designs to spice up your wardrobe.



Microfiber



Stylish



Patterns

Compression Level

20-30mmHg

Features and Benefits

- 100% of microfiber sits next to skin for a soft feel
- Fashionable, durable socks for everyday wear
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off
- Variety of patterns and colours available

Fibre Content

- 87% Nylon, 13% Spandex
(content can vary slightly by style)
- Latex free

Best Suited for

- Management of symptoms related to chronic venous disease
- Patients who prefer fashionable options
- Long-distance travel



Sizing Chart

	Small		Medium		Large	
Ankle Circumference	18–21.5cm (7"–8.5")		21.5–25.5cm (8.5"–10")		25.5–29cm (10"–11.5")	
Calf Circumference	28–39cm (11"–15.5")		31.5–45cm (12.5"–17.5")		35.5–51cm (14"–20")	
Calf Length	Short	Long	Short	Long	Short	Long
	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")
Size	SS	SL	MS	ML	LS	LL

Men's Shoe sizes: Short leg 7–11, Long leg 9.5–14 Women's Shoe Sizes: Short leg 5.5–9.5, Long leg 8–12

Styles and Colours

Closed Toe = ●								
Calf (C)	●	●	●	●	●	●	●	●
Closed Toe = ●								
Calf (C)	●	●	●	●	●	●		

Product Code Example: 832CMLW48

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour
83	2	C	ML	W	48

sigvaris

Motion Products

Previous	Series Number	New	Page
CUSHIONED COTTON	360	Cushioned Cotton	30
TRAVERSE SOCKS	412	High Tech	32
MERINO WOOL	422	Thermoregulating Wool	34
SWISS MOUNTAIN SOCK	462	Swiss Mountain Sock	36
DAILY COMFORT	850	Comfort	38

Motion

MOTION

Push Yourself

When daily comfort for hard-working legs is required, put your clients in motion. Sigvaris Motion socks offer **all day strong** support for people with active lives.

Why choose **Motion**?

- Cushioned sole for added comfort.
- Durable thermoregulating wool.
- Moisture-wicking yarns.

360

Motion CUSHIONED COTTON

for women & men

Motion Cushioned Cotton is soft and breathable with a padded sole for optimal comfort, and is ideal for sensitive skin.



**Cushioned
Sole**



Cotton



**Sensitive
Skin**

Compression Level

20-30mmHg

Features and Benefits

- Supima cotton for highest level of comfort and breathability
- Extra-padded sole, heel and toe to protect feet
- Flat, low-friction toe seam
- Diabetic friendly

Fibre Content

- 64% Supima Cotton, 28% Nylon, 8% Spandex
- Latex free

Best Suited for

- Management of CVD symptoms
- People who stand for prolonged periods of time
- Patients with sensitive skin and toes
- Work shoes and boots
- Casual and sports shoes







Sizing Chart

	Small		Medium		Large		X Large	
Ankle Circumference	18–21.5cm (7"–8.5")		21.5–25.5cm (8.5"–10")		25.5–29cm (10"–11.5")		29–36cm (11.5"–14")	
Calf Circumference	28–39cm (11"–15.5")		31.5–45cm (12.5"–17.5")		35.5–51cm (14"–20")		41–61cm (16"–24")	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")
Size	SS	SL	MS	ML	LS	LL	XS	XL

Women's Shoe Sizes: Short 5.5–9.5, Long 8–12 Men's Shoe Sizes: Short 7–11, Long 9.5–14

Styles and Colours

Closed Toe = ●		
		
	WHITE 00	BLACK 99
Calf (C)		

Product Code Example: 362CMLM00

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour
36	2	C	ML	M	00

Motion HIGH TECH

for women & men

Motion High Tech delivers optimal comfort where it counts with an anatomically-contoured shape and a cushioned sole.



**Moisture
Wicking**



**Anatomically
Contoured**



Colours

Compression Level

20-30mmHg

Features and Benefits

- High-tech fibres for improved thermoregulation
- Cushioned footbed, toe, and Achilles tendon zones for optimal comfort
- Additional compression in calf improves circulation

Fibre Content

- 64% Nylon, 17% Polyester, 19% Spandex
- Latex free

Best Suited for

- Management of symptoms related to chronic venous disease
- People who stand on their feet all day
- Walking and everyday activities
- Casual and occupational footwear
- Active individuals



Sizing Chart

	Small			Medium			Large			X Large		
Ankle Circumference	18–21.5cm (7"–8.5")			21.5–25.5cm (8.5"–10")			25.5–29cm (10"–11.5")			29–36cm (11.5"–14")		
Calf Circumference	28–39cm (11"–15.5")			31.5–45cm (11"–15.5")			35.5–51cm (14"–20")			40.5–58cm (16"–24")		
Women's Shoe Size	4.5–7.5	8–10.5	11–13.5	4.5–7.5	8–10.5	11–13.5	8–10.5	11–13.5	>14	8–10.5	11–13.5	>14
Men's Shoe Size	≤6.5	7–9.5	10–12.5	≤6.5	7–9.5	10–12.5	7–9.5	10–12.5	13–15.5	7–9.5	10–12.5	13–15.5
Size	SS	SM	SL	MS	MM	ML	LS	LM	LL	XS	XM	XL

Styles and Colours

Closed Toe = ●									
	WHITE 00	BLUE 50	LIME 54	PINK 56	RED 59	STEEL BLUE 62	BLOOD ORANGE 63	LIMEADE 64	BLACK 99
Calf (C)	●	●	●	●	●	●	●	●	●

Product Code Example: 412CML64

Series	mmHg Level	Style	Size	Colour
41	2	C	ML	64

Motion

THERMOREGULATING WOOL

for women & men

Motion Thermoregulating Wool combines the benefits of natural Merino wool with a cushioned sole, making it ideal for wherever life takes you!



Wool



Cushioned Sole



Thermoregulating

Compression Level

20-30mmHg

Features and Benefits

- Natural properties of wool make these socks odour-resistant
- Cushioned, padded soles help prevent blisters
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off

Fibre Content

- 67% Fine Merino Wool, 26% Nylon, 7% Spandex
- Latex free

Best Suited for

- Management of symptoms related to chronic venous disease
- People who stand on their feet all day
- Work boots, steel-toed shoes, and most occupational footwear
- Outdoor work environments



Sizing Chart

	Small			Medium			Large		
Ankle Circumference	18-21.5cm (7"-8.5")			21.5-25.5cm (8.5"-10")			25.5-29cm (10"-11.5")		
Calf Circumference	28-39cm (11"-15.5")			31.5-45cm (12.5"-17.5")			35.5-51cm (14"-20")		
Women's Shoe Size	4.5-7.5	8-10.5	11-13.5	4.5-7.5	8-10.5	11-13.5	8-10.5	11-13.5	>14
Men's Shoe Size	≤6.5	7-9.5	10-12.5	≤6.5	7-9.5	10-12.5	7-9.5	10-12.5	13-15.5
Size	SS	SM	SL	MS	MM	ML	LS	LM	LL

Styles and Colours

Closed Toe = ●	<div>CHARCOAL 12</div>	
	<div>OLIVE 31</div>	
Calf (C)	●	●

Product Code Example: 422CML12

Series	mmHg Level	Style	Size	Colour
42	2	C	ML	12

Motion

SWISS MOUNTAIN SOCK for women & men

The Swiss Mountain Sock is moisture-wicking and anatomically contoured, with a padded sole and Achilles tendon for ultimate comfort on the go.



**Cushioned
Sole**



**Anatomically
Contoured**



**Moisture
Wicking**

Compression Level

20-30mmHg

Features and Benefits

- Thermal and odour control fibres to promote moisture management
- Footbed, toe and Achilles tendon cushion zones help reduce blisters
- Swiss precision graduated compression

Fibre Content

- 52% Nylon, 25% Polypropylene, 12% Polyester, 11% Spandex
- Latex free

Best Suited for

- People who work on their feet all day
- Work boots, steel-toed shoes, and most occupational footwear
- Walking and outdoor activities



Sizing Chart

	Small			Medium			Large			X Large		
Ankle Circumference	18–21.5cm (7"–8.5")			21.5–25.5cm (8.5"–10")			25.5–29cm (10"–11.5")			29–33cm (11.5"–13")		
Calf Circumference	28–39cm (11"–15.5")			31.5–45cm (12.5"–17.5")			35.5–51cm (14"–20")			40.5–58cm (16"–24")		
Women's Shoe Size	4.5–7.5	8–10.5	11–13.5	4.5–7.5	8–10.5	11–13.5	8–10.5	11–13.5	>14	8–10.5	11–13.5	>14
Men's Shoe Size	≤6.5	7–9.5	10–12.5	≤6.5	7–9.5	10–12.5	7–9.5	10–12.5	13–15.5	7–9.5	10–12.5	13–15.5
Size	SS	SM	SL	MS	MM	ML	LS	LM	LL	XS	XM	XL

Women's Shoe Sizes: Short 5.5–9.5, Long 8–12 Men's Shoe Sizes: Short 7–11, Long 9.5–14

Styles and Colour

Closed Toe = ●



Calf (C)



Product Code Example: 462CSM59

Series	mmHg Level	Style	Size	Colour
46	2	C	SM	59

Motion COMFORT

for women & men

Motion Comfort offers cushioned soles and moisture wicking to keep your legs feeling energized all day long.



**Cushioned
Sole**



Durable



**Moisture
Wicking**

Compression Level

20-30mmHg

Features and Benefits

- NILIT® Aquarius fibres provide moisture wicking to keep feet dry and comfortable
- Cushioned soles to absorb impact and help prevent blisters
- Smooth, flat toe seam

Fibre Content

- 86% Nylon, 14% Spandex
- Latex free

Best Suited for

- Management of symptoms related to chronic venous disease
- Healthcare professionals
- Uniforms requiring black or solid colour
- Business travelers
- Individuals working in physically-demanding occupations
- Active individuals



Sizing Chart

	Small		Medium		Large	
Ankle Circumference	18–21.5cm (7"–8.5")		21.5–25.5cm (8.5"–10")		25.5–29cm (10"–11.5")	
Calf Circumference	28–39cm (11"–15.5")		31.5–45cm (12.5"–17.5")		35.5–51cm (10"–11.5")	
Calf Length	Short	Long	Short	Long	Short	Long
	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes: Short leg 5.5–9.5, Long leg 8–12 Men's Shoe sizes: Short leg 7–11, Long leg 9.5–14

Styles and Colour

Closed Toe = ●

BLACK

99

Calf (C)

●

Product Code Example: 852CMLM99


Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour
85	2	C	ML	M	99

sigvaris

Specialty Products

Product	Series Number	Page
Natural Rubber	500	42
Advance Armsleeves	910	44

Specialty



SPECIALTY

Some patients may require a more specialized product to help relieve the effects of more advanced edemas and other complications of venous disorders.

SIGVARIS GROUP offers both upper and lower extremity specialty solutions that have a higher level of containment and durability while providing the same comfort.

Why choose **Specialty**?

- Specific medical needs.
- Need low fatigue fibers.
- Natural rubber.

500

Specialty NATURAL RUBBER

for women & men

Natural Rubber offers patients a durable, long-term high-containment solution for compression therapy with minimal fibre fatigue throughout the day.



Natural Rubber



Durable



Wide Size Range

Compression Level

30-40mmHg

40-50mmHg

50-60mmHg

Features and Benefits

- Ideal for patients with more serious indications
- Unique sizing for a custom-like fit
- Comfortable grip-top band on thigh styles to prevent slipping
- Caution: contains natural rubber latex that may cause allergic reaction

Fibre Content

- 70% Nylon, 30% Natural Latex Rubber
- 75% Nylon, 20% Natural Latex Rubber, 5% other fibres (Pantyhose style)

Best Suited for

- Severe varicosities
- Chronic edema post DVT
- Severe venous stasis
- Prevention of ulcer recurrence
- Early stages of lymphedema
- Active venous ulcers



Silicone
dotted grip-top
band in thigh

Sizing Chart

	Small				Medium				Large				X Large			
Ankle Circumference	19–23cm (7.5"–9")				23–30cm (9"–12")				30–40cm (12"–15.5")				40–46cm (15.5"–18")			
Calf Circumference	Average		Full		Average		Full		Average		Full		Average		Full	
	28–38cm (11"–15")		<43cm (17")		32–43cm (12.5"–17")		<48cm (19")		36–48cm (14.5"–19")		<57cm (22.5")		45–57cm (18"–22.5")		68cm (27")	
Calf Circumference for 505C (Full only)	505C: full calf only		32–43cm (12.5"–17")		505C: full calf only		36–48cm (14.5"–19")		505C: full calf only		45–57cm (18"–22.5")		505C: not available in XL			
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long
	40cm (≤15.5")	40cm (≤15.5")	40cm (≤15.5")	40cm (≤15.5")	40cm (≤15.5")	40cm (≤15.5")	40cm (≤15.5")	40cm (≤15.5")	40cm (≤15.5")	40cm (≤15.5")	40cm (≤15.5")	40cm (≤15.5")	40cm (≤15.5")	40cm (≤15.5")	40cm (≤15.5")	40cm (≤15.5")
Thigh Circumference	Average		Full		Average		Full		Average		Full		N/A			
	41–54cm (16"–21")		<61cm (24")		48–61cm (19"–24")		≤71cm (28")		58–71cm (23"–28")		≤81cm (32")					
Thigh Circumference for 505T (Full only)	505T: full calf only		<61cm (24")		505T: full calf only		≤71cm (28")		505T: full calf only		≤81cm (32")					
Leg Length	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long				
	≤71cm (28")	>71cm (28")	≤71cm (28")	>71cm (28")	≤71cm (28")	>71cm (28")	≤71cm (28")	>71cm (28")	≤71cm (28")	>71cm (28")	≤71cm (28")	>71cm (28")				
Hip Circumference Guidelines	<122cm (48")				<137cm (48")				<155cm (61")							
Size	S1	S2	S3	S4	M1	M2	M3	M4	L1	L2	L3	L4	X1	X2	X3	X4

30–40mmHg: all sizes except XL available

40–50mmHg: only available in full sizes

Styles and Colour

Open Toe = O



Caution: this product contains natural rubber latex which may cause allergic reactions. For all indications & contraindications, please see Contraindications & Cautions on page 87. Sigvaris accepts no liability for non-observance of contraindications or cautions.

Calf (C)	O	Calf only available in S, M, and L full circumference for 50–60mmHg
Thigh w/out Grip-Top (T)	O	
Thigh w/Grip-Top (N) Silicone grip-top band	O	Not available in 50–60mmHg
Pantyhose (P)	O	Pantyhose only available in 30–40mmHg and only in full thigh circumference
Thigh w/Waist Attachment (W) Indicate leg LEFT (/L) or RIGHT (/R)	O	Thigh with waist attachment only available in 30–40mmHg and 40–50mmHg and only available in full thigh circumference

Product Code Example: 504CS2O77

Series	mmHg Level	Style	Size	Gender/Open Toe	Colour
50	4	C	S2	O	77

Product Code Example: 503WL1O77/L

Series	mmHg Level	Style	Size	Gender/Open Toe	Colour	Leg
50	3	W	L1	O	77	/L

Specialty ADVANCE ARMSLEEVE

for women & men

The Advance Armsleeve combines the ultimate softness with silver ion properties and a Sensinnov grip-top band that's ideal for all-day comfort and effectiveness.



**Moisture
Wicking**



**Sensinnov
Grip-top
Band**



Breathable

Compression Level

15-20mmHg

20-30mmHg

30-40mmHg

Features and Benefits

- Ideal for patients with arm lymphedema or post-mastectomy care
- Sensinnov top-band ensures armsleeve stays in place all day long
- Provides moisture wicking properties to keep skin dry and comfortable
- Breathable to provide comfortable, year-round wear

Fibre Content

- 50% Nylon, 28% Rayon, 22% Spandex
(Content can vary slightly by style)
- Latex free

Best Suited for

- Early stages of lymphedema
- Patients with arm lymphedema
- Post-mastectomy care



**Sensinnov
grip-top
band**



Sizing Chart

	Small		Medium				Large				X Large			
Palm	17–21cm (6.5"–8")		18–22cm (7"–8.5")				19–23cm (7.5"–9")				20–24cm (8"–9.5")			
Wrist	14–16cm (5.5"–6.5")		16–18cm (6.5"–7")				18–20cm (7"–8")				20–23cm (8"–9")			
Upper Mid-Arm	Regular		Regular		Plus		Regular		Plus		Regular		Plus	
	26–34cm (10"–13.5")		29–37cm (11.5"–14.5")		37–43cm (14.5"–17")		32–40cm (12.5"–15.5")		40–46cm (15.5"–18")		34–42cm (13.5"–16.5")		42–48cm (16.5"–19")	
Arm Length	Regular	Long	Regular	Long	Regular	Long	Regular	Long	Regular	Long	Regular	Long	Regular	Long
	≤43cm (≤17")	>43cm (>17")	≤43cm (≤17")	>43cm (>17")	≤43cm (≤17")	>43cm (>17")	≤43cm (≤17")	>43cm (>17")	≤43cm (≤17")	>43cm (>17")	≤43cm (≤17")	>43cm (>17")	≤43cm (≤17")	>43cm (>17")
Size	SR	SL	MR	ML	MR/P	ML/P	LR	LL	LR/P	LL/P	XR	XL	XR/P	XL/P

Styles and Colour

Sleeve = ●




Armsleeve
w/Gauntlet (G)

Sensinnov grip-top band



Armsleeve
w/out Gauntlet (O)

Sensinnov grip-top band



Note: all armsleeves come with Grip-top (/S)

15–20mmHg not available with gauntlet
and not available in plus sizing

Armsleeve without gauntlet

Product Code Example: 912ASLO77/S

Series	mmHg Level	Style	Size	Gauntlet	Colour	Grip-Top
91	2	A	SL	O	77	/S

Armsleeve with gauntlet

Product Code Example: 912ALRG77/PS

Series	mmHg Level	Style	Size	Gauntlet	Colour	Plus Size	Grip-Top
91	2	A	LR	G	77	/P	S

CONSUMER HEALTH CARE PRODUCTS

Product	Series Number	Page
CASUAL COTTON	146 / 186	48
CUSHIONED COTTON	142 / 182	50
SEA ISLAND COTTON	151 / 191	52
ALL-SEASON MERINO WOOL	152 / 192	54
MICROFIBER SHADES	143 / 183	56
SHEER FASHION	120	58
SOFT SILHOUETTE LEGGINGS	170	60
BUSINESS CASUAL	189	62
ATHLETIC RECOVERY SOCKS	401	64
MERINO OUTDOOR SOCKS	421	66
EVERSOFT DIABETIC SOCKS	160	68
PERFORMANCE SLEEVES	412	70
TRAVENO TRAVEL SOCKS	146/T & 186/T	72



CONSUMER HEALTH CARE

Improve leg health with a simple change of socks.

For many people, experiencing tired or achy legs is the first step in learning more about the health benefits of graduated compression.

The SIGVARIS CONSUMER HEALTH CARE line introduces many of the therapeutic benefits of graduated compression with stylish options in a simplified sizing system. Consumers will be able to find their new favourite pair of socks or hosiery with SIGVARIS.

SIGVARIS

146/186

CASUAL COTTON

for women & men

Compression Level

15-20mmHg

Features and Benefits

- » Made with premium Supima Cotton for optimum softness and durability
- » Flat, low-friction toe seam to prevent irritation
- » Soft, stretchable knit-in band prevents pinching

Fibre Content

- » 66% Supima Cotton, 26% Nylon, 8% Spandex
- » Latex free

Best Suited for

- » Tired, heavy, restless, aching legs
- » Travelers, occupational leg health risks, expecting mothers, family history of chronic venous disorder



Sizing Chart				
146	Women's Shoe Size	5-7	7.5-9.5	10-12
	Ankle Circumference	18-23cm (7"-9")	19-25cm (7.5"-10")	20-27cm (8"-10.5")
186	Men's Shoe Size	6.5-8.5	9-11	11.5-14
	Ankle Circumference	19-25cm (7.5"-10")	20-27cm (7.5"-10")	21-29cm (8.5"-11.5")
Size		A	B	C

Style and Colours					
Closed Toe = ●					
	WHITE 00	NAVY 10	BROWN 11	KHAKI 30	BLACK 99
Calf (C)					

Product Code Example: **146CB10**

Series	Style	Size	Colour
146	C	B	10



146
Women's Calf
15-20mmHg
Calf Size B in Khaki

146CB30



186
Men's Calf
15-20mmHg
Calf Size C in Black

186CC99

CUSHIONED COTTON

for women & men

Compression Level

15–20mmHg

Also available in 20–30mmHg, see page 30

Features and Benefits

- » Cushioned sole, heel and toe area to protect feet
- » Flat, low-friction toe seam
- » Soft, stretchable knit-in calf band prevents pinching
- » Ideal for everyday wear and for sensitive skin
- » Suitable for patients with diabetes
- » No compression in heel and instep for ease of donning

Fibre Content

- » 68% Supima Cotton, 25% Nylon, 7% Spandex
- » Latex-Free

Best Suited for

- » Tired, heavy, restless, aching legs
- » Travelers, occupational leg health risks, expecting mothers, family history of chronic venous disorders



Sizing Chart				
142	Women's Shoe Size	5-7	7.5-9.5	10-12
	Ankle Circumference	18-23cm (7"-9")	19-25cm (7.5"-10")	20-27cm (8"-10.5")
182	Men's Shoe Size	6.5-8.5	9-11	11.5-14
	Ankle Circumference	19-25cm (7.5"-10")	20-27cm (8"-10.5")	21-29cm (8.5"-11.5")
Size		A	B	C

Style and Colours

Closed Toe = ●

WHITE
00

BLACK
99

Calf (C)

●

●

Product Code Example: **142CB00**

Series	Style	Size	Colour
142	C	B	00



15-20mmHg
Calf Size B in White

142CB00



15-20mmHg
Calf Size C in Black

182CC99

SEA ISLAND COTTON

for women & men

Compression Level

15-20mmHg

Also available in 20-30mmHg, see page 20

Features and Benefits

- » Knitted from the world's rarest cotton
- » Sea Island Cotton yarns are lightweight and have great strength, a silky luster, and an incredibly soft feel
- » Breathable and comfortable for sensitive skin
- » Ideal for travel and everyday wear

Fibre Content

- » 43% Sea Island Cotton, 39% Nylon, 11% Spandex, 7% Supima Cotton
- » Latex-Free

Best Suited for

- » Tired, heavy, restless, aching legs
- » Travelers, occupational leg health risks, expecting mothers, family history of chronic venous disorder



ZURICH
COLLECTION

Sizing Chart				
151	Women's Shoe Size	5-7	7.5-9.5	10-12
	Ankle Circumference	18-23cm (7"-9")	19-25cm (7.5"-10")	20-27cm (8"-10.5")
191	Men's Shoe Size	6.5-8.5	9-11	11.5-14
	Ankle Circumference	19-25cm (7.5"-10")	20-27cm (8"-10.5")	21-29cm (8.5"-11.5")
Size		A	B	C

Style and Colours			
Closed Toe = ●			
<div>NAVY 10</div> <div>BROWN 11</div> <div>BLACK 99</div>			
Calf (C)	●	●	●

Product Code Example: 151CA11

Series	Style	Size	Colour
151	C	A	11



15-20mmHg
Calf Size A in Navy

151CA10



15-20mmHg
Calf Size B in Black

191CB99

ALL SEASON MERINO WOOL

for women & men

Compression Level

15-20mmHg

Also available in 20-30mmHg, see page 22

Features and Benefits

- » Soft, comfortable and itch-free
- » Natural properties of wool make these socks odour resistant
- » Thermoregulating for year-round wear and keeps feet dry

Fibre Content

- » 53% Extra-fine Merino Wool, 38% Nylon, 9% Spandex
- » Latex-Free

Best Suited for

- » Tired, heavy, restless, aching legs
- » Travelers, occupational leg health risks, expecting mothers, family history of chronic venous disorder



ZURICH
COLLECTION

Sizing Chart				
152	Women's Shoe Size	5-7	7.5-9.5	10-12
	Ankle Circumference	18-23cm (7"-9")	19-25cm (7.5"-10")	20-27cm (8"-10.5")
192	Men's Shoe Size	6.5-8.5	9-11	11.5-14
	Ankle Circumference	19-25cm (7.5"-10")	20-27cm (8"-10.5")	21-29cm (8.5"-11.5")
Size		A	B	C

Style and Colours			
Closed Toe = ●		<div>NAVY 10</div>	<div>BROWN 11</div>
			<div>BLACK 99</div>
Calf (C)	●	●	●

Product Code Example: 152CA11

Series	Style	Size	Colour
152	C	A	11



15-20mmHg
Calf Size A in Black
152CA99



15-20mmHg
Calf Size C in Navy
192CC10

MICROFIBER SHADES

for women & men

Compression Level

15-20mmHg

Also available in 20-30mmHg, see page 26

Features and Benefits

- » Advanced microfiber technology for ultimate comfort
- » Glides on and off with ease
- » Excellent breathability and moisture management
- » Soft and durable sock for everyday wear
- » Variety of patterns and colours available

Fibre Content

- » 87% Nylon, 13% Spandex
- » Latex-Free

Best Suited for

- » Tired, heavy, restless, aching legs
- » Travelers, occupational leg health risks, expecting mothers, family history of chronic venous disorder



Sizing Chart				
143	Women's Shoe Size	5-7	7.5-9.5	10-12
	Ankle Circumference	18-23cm (7"-9")	19-25cm (7.5"-10")	20-27cm (8"-10.5")
183	Men's Shoe Size	6.5-8.5	9-11	11.5-14
	Ankle Circumference	19-25cm (7.5"-10")	20-27cm (8"-10.5")	21-29cm (8.5"-11.5")
Size		A	B	C

Style and Colours								
Closed Toe = ● <i>Actual colour may vary</i>								
								
Calf (C)		●	●	●	●	●	●	●
Closed Toe = ● <i>Actual colour may vary</i>								
Calf (C)		●	●	●	●	●	●	

Product Code Example: **143CA46**

Series	Style	Size	Colour
143	C	A	46



15-20mmHg
Calf Size B in Purple Argyle

143CB46



15-20mmHg
Calf Size C in Graphite Mini-Stripe

183CC07

SHEER FASHION *for women*

Compression Level

15-20mmHg

Features and Benefits

- » Fabric feels soft against the skin
- » Incredible sheer design in a variety of fashionable colours
- » Sheer toe and heel area
- » Comfortable stretch is excellent for everyday wear
- » Cotton-lined gusset on pantyhose styles

Fibre Content

- » 79% Nylon, 21% Spandex
- » Content can vary slightly by style
- » Latex-Free

Best Suited for

- » Tired, heavy, restless, aching legs
- » Travelers, occupational leg health risks, expecting mothers, family history of chronic venous disorder



15-20mmHg
Closed Toe Calf Size B in Taupe

120CB29



15-20mmHg
Closed Toe Pantyhose Size C in Black

120PC99

Pantyhose and Maternity Sizing Chart

Cm/Ft	Kg/Lbs	45 / 100	48 / 105	50 / 110	52 / 115	54 / 120	57 / 125	59 / 130	61 / 135	64 / 140	66 / 145	68 / 150	70 / 155	73 / 160	75 / 165	77 / 170	80 / 175	82 / 180	84 / 185	86 / 190	89 / 195	91 / 200	93 / 205	95 / 210	98 / 215	100 / 220	102 / 225	105 / 230
145cm / 4'9"																												
147cm / 4'10"																												
150cm / 4'11"																												
152cm / 5'0"																												
155cm / 5'1"																												
157cm / 5'2"																												
160cm / 5'3"																												
162cm / 5'4"																												
165cm / 5'5"																												
167cm / 5'6"																												
170cm / 5'7"																												
173cm / 5'8"																												
175cm / 5'9"																												
178cm / 5'10"																												
180cm / 5'11"																												
183cm / 6'0"																												

Note: please use pre-pregnancy weight for maternity pantyhose. * Maternity pantyhose not available in size F.

Calf / Thigh Sizing Chart

Women's Shoe Size	5-7	7.5-9.5	10-12
Ankle Circumference	18-23cm (7"-9")	19-25cm (7.5"-10")	20-27cm (8"-10.5")
Size	A	B	C

Weight Guidelines For Thigh Style: Size A: <135Lbs., Size B: <180Lbs., Size C: <230Lbs.

Styles and Colours

							
	CHARCOAL 12	TAUPE 29	NATURAL 33	SUNTAN 36	MULBERRY 60	CAFÉ 73	BLACK 99
Calf (C)	●	●●	●●	●●	●	●	●
Thigh (N)	●	●	●	●		●	●
Pantyhose (P)	●	●	●	●		●	●
Maternity Pantyhose (M)			●				●

Product Code Example: 120PE73

Series	Style	Size	Colour
120	P	E	73

Product Code Example: 120CCO36

Series	Style	Size	Open Toe	Colour
120	C	C	O	36

SOFT SILHOUETTE LEGGINGS

for women

Compression Level

15-20mmHg

Features and Benefits

- » Stylish, subtle ribbed design looks chic and pairs perfectly with dresses, blouses, tunics and shirts
- » Comfort-stretch waist band provides a seamless look for all-day wearing comfort
- » Advanced soft stretch microfiber fabric shapes legs and smooths skin

Fibre Content

- » 83% Nylon, 17% Spandex
- » Latex-Free

Best Suited for

- » Tired, heavy, restless, aching legs,
- » Occupational leg health risks, family history of chronic venous disorders



Sizing Chart															
cm / ft	Kg / Lbs														
	45 / 100	48 / 105	50 / 110	52 / 115	54 / 120	57 / 125	59 / 130	61 / 135	64 / 140	66 / 145	68 / 150	70 / 155	73 / 160	75 / 165	77 / 170
147cm / 4'10"															
150cm / 4'11"															
152cm / 5'0"															
155cm / 5'1"															
157cm / 5'2"															
160cm / 5'3"															
162cm / 5'4"															
165cm / 5'5"															
167cm / 5'6"															
170cm / 5'7"															
173cm / 5'8"															
175cm / 5'9"															
178cm / 5'10"															
180cm / 5'11"															
183cm / 6'0"															

Guidelines for Ankle Circumference
Size A: 7"-9", Size B & C: 8"-10", Size D: 8.5"-11.5"

Style and Colours

Open Foot = ●

BLACK 99

Leggings (L)

●

Product Code Example: 170LB99			
Series	Style	Size	Colour
170	L	B	99



15-20mmHg
Legging Size B in Black

170LB99



Non-binding ankle cuff
stays comfortably in place

BUSINESS CASUAL

for men

Compression Level

15-20mmHg

Features and Benefits

- » Stretchable toe fabric and heel design increase comfort and durability
- » Classic ribbed look for everyday wear
- » Wide, non-binding top band for comfortable fit
- » Soft, microfiber inlay yarns

Fibre Content

- » 88% Nylon, 12% Spandex
- » Latex-Free

Best Suited for

- » Tired, heavy, restless, aching legs,
- » Travelers, occupational leg health risks, family history of chronic venous disorders



Sizing Chart

Men's Shoe Size	6.5–8.5	9–11	11.5–14
Ankle Circumference	19–25cm (7.5"–10")	20–27cm (8"–10.5")	21–29cm (8.5"–11.5")
Size	A	B	C

Style and Colours

Closed Toe = ●	<div>NAVY 10</div>	<div>BROWN 11</div>	<div>CHARCOAL 12</div>	<div>BLACK 99</div>
Calf (C)	●	●	●	●

Product Code Example: 189CB12

Series	Style	Size	Colour
189	C	B	12



15–20mmHg
Calf Size B in Charcoal

189CB12



15–20mmHg
Calf Size C in Brown

189CC11

ATHLETIC RECOVERY SOCKS

for women & men

Compression Level

15-20mmHg

Features and Benefits

- » Helps reduce lactic acid buildup
- » Comfortable and easy to put on and take off
- » *drirelease* yarns provides moisture wicking and odour control
- » No compression in heel and instep for ease of donning

Fibre Content

- » 67% *drirelease* Polyester (57% Polyester, 10% Cotton), 26% Nylon, 7% Spandex
- » Latex-Free

Best Suited for

- » Tired, heavy, restless, aching legs
- » Post-workout recovery
- » Mild athletic activity, walking



Sizing Chart					
Women's Shoe Size	<7.5	7.5–9.5	10–12	12.5–15	≥15.5
Men's Shoe Size	<6.5	6.5–8.5	9–11	11.5–14	14.5–17
Ankle Circumference	18–24cm (7"–9.5")	19–26cm (7.5"–10")	20–27cm (8"–10.5")	21–29cm (8.5"–11.5")	21–29cm (8.5"–11.5")
Size	S	M	L	X	2

Style and Colours

Closed Toe = ●

WHITE
00

BLACK
99

Calf (C)

●

●

Product Code Example: **401C200**

Series	Style	Size	Colour
401	C	2	00



15–20mmHg
Calf Size XX-Large in Black

401C299



15–20mmHg
Calf Size Medium in White

401CM00

MERINO OUTDOOR SOCKS

for women & men

Compression Level

15–20mmHg

Also available in 20–30mmHg, see page 34

Features and Benefits

- » Soft, comfortable and itch-free
- » Natural properties of wool make these socks odour-resistant
- » Mesh flex zones for comfort and ventilation
- » Cushioned, padded soles protect against blisters
- » No compression in heel and instep for ease of donning

Fibre Content

- » 67% Fine Merino Wool, 26% Nylon, 7% Spandex
- » Latex-Free

Best Suited for

- » Tired, heavy, restless, aching legs
- » Hiking, fishing & hunting
- » Outdoor sports and walking



Sizing Chart				
Women's Shoe Size	<7.5	7.5–9.5	10–12	12.5–15
Men's Shoe Size	<6.5	6.5–8.5	9–11	11.5–14
Ankle Circumference	18–24cm (7"–9.5")	19–26cm (7.5"–10")	20–27cm (8"–10.5")	21–29cm (8.5"–11.5")
Size	S	M	L	X

Style and Colours		
Closed Toe = ●	CHARCOAL 12	OLIVE 31
	●	●
Calf (C)		

Product Code Example: **421CS31**

Series	Style	Size	Colour
421	C	S	31



15–20mmHg
Calf Size Large in Olive

421CL31



15–20mmHg
Calf Size Medium in Charcoal

421CM12

EVERSOFT DIABETIC SOCKS

for women & men

Compression Level

8-15mmHg

Features and Benefits

- » Thick, padded soles, heels, and toes for comfort and protection
- » Flat, low friction toe seam
- » Soft, stretchable knit-in band prevents pinching
- » Compression from ankle to knee, compression-free foot

Fibre Content

- » 67% *drirelease* Polyester (57% Polyester, 10% Cotton), 26% Nylon, 7% Spandex
- » Latex-Free

Best Suited for

- » Individuals with diabetes
- » Individuals with limited dexterity or lower tolerance for compression
- » Individuals with sensitive skin



Sizing Chart				
Women's Shoe Size	<7.5	7.5–9.5	10–12	12.5–15
Men's Shoe Size	<6.5	6.5–8.5	9–11	11.5–14
Ankle Circumference	18–24cm (7"–9.5")	19–26cm (7.5"–10")	20–27cm (8"–10.5")	21–29cm (8.5"–11.5")
Size	S	M	L	X

Style and Colours

Closed Toe = ●



Calf (C)

Product Code Example: 160CS00

Series	Style	Size	Colour
160	C	S	00



8–15mmHg
Calf Size Small in White

160CS00



8–15mmHg
Calf Size Large in White

160CL00

PERFORMANCE SLEEVES

for women & men

Compression Level

20-30mmHg

Features and Benefits

- » Fabric resists deterioration in salt, chlorinated or fresh water
- » Stabilizes and diminishes vibrations of active muscles
- » Quick drying
- » High-visibility reflective transfer

Fibre Content

- » 75% Nylon, 25% Spandex
- » Latex-Free

Best Suited for

- » Running, triathlons, endurance sports and gym workouts
- » Team sports



Sizing Chart

Ankle Circumference	18–21.5cm (7"–8.5")	21.5–25.5cm (8.5"–10")	25.5–29cm (10"–11.5")	29–33cm (11.5"–13")
Calf Circumference	28–39cm (11"–15.5")	31.5–45cm (12.5"–17.5")	35.5–51cm (14"–20")	40.5–58cm (16"–23")
Size	S	M	L	X

Style and Colours

Calf Sleeve = ●	WHITE 00	BLUE 50	LIME 54	PINK 56*	PURPLE 58*	RED 59	STEEL BLUE 62	BLOOD ORANGE 63	LIMEADE 64	BLACK 99
Calf Sleeve (V)	●	●	●	●	●	●	●	●	●	●

* Pink and purple not available in X-Large

Product Code Example: 412VL63

Series	Style	Size	Colour
412	V	L	63



20–30mmHg
Calf Sleeve Size Large in Blood Orange

412VL63



20–30mmHg
Calf Sleeve Size Medium in Steel Blue

412VM62

NEW!

TRAVENO TRAVEL SOCKS

for women & men

Compression Level

15–20mmHg

Features and Benefits

- » Made with premium Supima Cotton for optimum softness and durability
- » Flat, low-friction toe seam to prevent irritation
- » Soft, stretchable knit-in band prevents pinching

Fibre Content

- » 66% Supima Cotton, 26% Nylon, 8% Spandex
- » Latex free

Best Suited for

- » Travelers
- » Tired, heavy, restless, aching legs
- » Occupational leg health risks, expecting mothers, family history of chronic venous disorders



Style and Colours

Closed Toe = ●

BLACK
99

Calf (C)



Sizing Chart

146/T	Women's Shoe Size	5–7	7.5–9.5	10–12
	Ankle Circumference	18–23cm (7"–9")	19–25cm (7.5"–10")	20–27cm (8"–10.5")
186/T	Men's Shoe Size	6.5–8.5	9–11	11.5–14
	Ankle Circumference	19–25cm (7.5"–10")	20–27cm (8"–10.5")	21–29cm (8.5"–11.5")
Size		A	B	C

IMPROVE YOUR
CIRCULATION
IN THE AIR
OR ON
THE GROUND

Traveno®
TRAVEL SOCKS



Travel Socks for Women & Men

Available in 15–20mmHg in black

Countertop displays available. Contact your Territory Manager or call **1.800.363.4999** for more information

sigvaris

Inelastic Wraps Products

Product	Page
Compreflex Standard Calf & Foot	76
Compreflex Standard Calf	76
Compreflex Transition Calf	78
Compreflex Complete Calf	78
Assorted Liners	80
Accutab	81
Compreflex Strap Extender	81
Assorted Compreboots	82

Inelastic Wraps



INELASTIC WRAPS

The Sigvaris line of products for medically complex extremities (MCE) is for the treatment and management of conditions impacting the lymphatic and venous systems.

The product portfolio contains ready-to-wear garments for the active treatment of lymphedema and for patients managing chronic venous diseases in the upper and lower extremities.

SIGVARIS is looking towards a brighter future with new therapies for complex venous and lymphatic diseases. Our inelastic wraps deliver maximum therapeutic results with increased patient compliance.

Why choose **Inelastic Wraps**?

- Adjustable pressure.
- Large variety of solutions.

Compreflex Standard **CALF & FOOT**

20-50mmHg



When to Use

- Adjustable ranges of resting compression is desired (20-50mmHg)
- Non-attached foot piece is desired
- No compression in the liner is desired
- Patient requires a calf and foot garment

Conditions

- Lymphedema
- Edema
- Chronic venous disease (CVD)
- Chronic venous insufficiency (CVI) including, but not limited to ulcer due to venous stasis
- Lipedema



**Compreboot
Standard**

30-40mmHg
(see page 82 for details)



Basic Liner

Non-Compressive
(see page 80 for details)



Accutab™ Devices

Accutab Devices
(see page 81 for details)

Compreflex Standard **CALF**

20-50mmHg



When to Use

- Varying ranges of adjustable resting compression is desired (20-50mmHg)
- No compression in the liner is desired
- No Compreboot needed
- Patient requires compression in the calf only or a different foot option (size or boot)
- Ability to wear standard footwear over liner is desired


Conditions

- Lymphedema
- Edema
- Chronic venous disease (CVD)
- Chronic venous insufficiency (CVI) including, but not limited to ulcer due to venous stasis
- Lipedema



Basic Liner

Non-Compressive
(see page 80 for details)



Accutab Devices


Accutab Devices
(see page 81 for details)

COMPREFLEX WRAPS



30–40mmHg Compreboot for added foot and ankle compression

REQUIRED MEASUREMENTS



Available Styles and Colours

	Black	Beige	Heathered Silver
	(99)	(77)	(N/A)
Compreflex Standard – Calf & Foot	●	●	
Compreflex Standard – Calf	●	●	
Basic Liner			●

Compreflex Standard Calf & Foot Sizing Chart (A–C1)

	Small	Medium	Large	X-Large	XX-Large
C1	29–39cm	34–44cm	40–50cm	46–56cm	54–64cm
C	29–39cm	34–44cm	40–50cm	46–56cm	54–64cm
B	24–34cm	29–39cm	34–44cm	39–49cm	44–55cm
A	19–26cm	21–30cm	26–36cm	31–41cm	36–46cm

	Black					
		Regular	Tall	Regular	Tall	Regular
Beige	Black	1401-BKR	1402-BKR	1403-BKR	1404-BKR	1405-BKR
		1401-BKT	1402-BKT	1403-BKT	1404-BKT	1405-BKT
Beige	Beige	1411-BKR	1412-BKR	1413-BKR	1414-BKR	1415-BKR
		1411-BKT	1412-BKT	1413-BKT	1414-BKT	1415-BKT

Compreflex Standard – Calf Sizing Chart (A–C1)

	Small	Medium	Large	X-Large	XX-Large
C1	29–39cm	34–44cm	40–50cm	46–56cm	54–64cm
C	29–39cm	34–44cm	40–50cm	46–56cm	54–64cm
B	24–34cm	29–39cm	34–44cm	39–49cm	44–55cm
A	19–26cm	21–30cm	26–36cm	31–41cm	36–46cm

	Black					
		Regular	Tall	Regular	Tall	Regular
Beige	Black	1401-NF-BKR	1402-NF-BKR	1403-NF-BKR	1404-NF-BKR	1405-NF-BKR
		1401-NF-BKT	1402-NF-BKT	1403-NF-BKT	1404-NF-BKT	1405-NF-BKT
Beige	Beige	1411-NF-BKR	1412-NF-BKR	1413-NF-BKR	1414-NF-BKR	1415-NF-BKR
		1411-NF-BKT	1412-NF-BKT	1413-NF-BKT	1414-NF-BKT	1415-NF-BKT

Compreflex Transition **CALF**

20-50mmHg



When to Use

- Patient requires a compressive liner; the Transition Liner provides 10-15mmHg in the foot and ankle
- Mild foot compression in liner is desired
- Ability to wear standard footwear over liner is desired

Conditions

- Lymphedema
- Edema
- Chronic venous disease (CVD)
- Chronic venous insufficiency (CVI) including, but not limited to ulcer due to venous stasis
- Lipedema



Transition Liner

10-15 mmHg of mild compression in the foot and ankle.
(see page 80 for details)



Accutab Devices

Accutab Devices
(see page 81 for details)



Compreflex Complete **CALF**

20-50mmHg



When to Use

- Patient requires a compressive liner; the Complete Liner provides 10-15mmHg in the foot and ankle
- Mild foot and calf compression in liner is desired
- Ability to wear standard footwear over liner is desired

Conditions

- Lymphedema
- Edema
- Chronic venous disease (CVD)
- Chronic venous insufficiency (CVI) including, but not limited to ulcer due to venous stasis
- Lipedema



Complete Liner

Mild compression in foot, ankle, and calf
(see page 80 for details)




Accutab Devices

Accutab Devices
(see page 81 for details)



COMPREFLEX WRAPS

REQUIRED MEASUREMENTS



Compreflex Complete Calf
20–50mmHg
Medium Regular in Beige
1412-FC-BKR

Available Styles and Colours

	Black	Beige	Heathered Silver
	(99)	(77)	(N/A)
Compreflex Transition – Calf	●	○	
Compreflex Complete – Calf	●	○	
Basic Liner			○

Liner that ships with wrap will match wrap colour

Compreflex Transition – Calf Sizing Chart (A–C1)

	Small	Medium	Large	X-Large	XX-Large
C1	29–39cm	34–44cm	40–50cm	46–56cm	54–64cm
C	29–39cm	34–44cm	40–50cm	46–56cm	54–64cm
B	24–34cm	29–39cm	34–44cm	39–49cm	44–55cm
A	19–26cm	21–30cm	26–36cm	31–41cm	36–46cm

Black	Regular	1401-UC-BKR	1402-UC-BKR	1403-UC-BKR	1404-UC-BKR	1405-UC-BKR
	Tall	1401-UC-BKT	1402-UC-BKT	1403-UC-BKT	1404-UC-BKT	1405-UC-BKT
Beige	Regular	1411-UC-BKR	1412-UC-BKR	1413-UC-BKR	1414-UC-BKR	1415-UC-BKR
	Tall	1411-UC-BKT	1412-UC-BKT	1413-UC-BKT	1414-UC-BKT	1415-UC-BKT

Compreflex Complete – Calf Sizing Chart (A–C1)

	Small	Medium	Large	X-Large	XX-Large
C1	29–39cm	34–44cm	40–50cm	46–56cm	54–64cm
C	29–39cm	34–44cm	40–50cm	46–56cm	54–64cm
B	24–34cm	29–39cm	34–44cm	39–49cm	44–55cm
A	19–26cm	21–30cm	26–36cm	31–41cm	36–46cm

Black	Regular	1401-FC-BKR	1402-FC-BKR	1403-FC-BKR	1404-FC-BKR	1405-FC-BKR
	Tall	1401-FC-BKT	1402-FC-BKT	1403-FC-BKT	1404-FC-BKT	1405-FC-BKT
Beige	Regular	1411-FC-BKR	1412-FC-BKR	1413-FC-BKR	1414-FC-BKR	1415-FC-BKR
	Tall	1411-FC-BKT	1412-FC-BKT	1413-FC-BKT	1414-FC-BKT	1415-FC-BKT

Basic Liner

Non-compressive

Features

- Non-compressive liner for the protection of the skin
- Easy to apply knee-high liners
- *drirelease* moisture wicking and odour control
- Flat seams for added comfort
- Silver ionized treatment for bacteriostatic properties
- Grey colour only

Size: One size fits all	
Pair	10 Pair Pack
1700 - NK1	1700 - NKC
Included With	
Compreflex Standard Calf & Foot Compreflex Standard Calf	



Transition Liner

10-15mmHg

Features

- Mild compression in foot and ankle only
- *drirelease* moisture wicking and odour control
- Extra thick sole from toe to heel for added comfort
- Silver ionized treatment for bacteriostatic properties
- Available in Black and Beige

Transition Liner Sizing Chart (Included with Compreflex Transition)

SIZE	Ankle Circ. (A)	Calf Circ. (C)	Women's Shoe	Men's Shoe	Pair (Black)	Pair (Beige)
Small	19-26cm	29-39cm	<7½	<6½	1601-CK	1611-CK
Medium	21-30cm	34-44cm	7.5-9.5	6.5-8.5	1602-CK	1612-CK
Large	26-36cm	40-50cm	10-12	9-11	1603-CK	1613-CK
X-Large	31-41cm	46-56cm	12.5-15	11.5-14	1604-CK	1614-CK
XX-Large	36-46cm	54-64cm	12.5-15	11.5-14	1605-CK	1615-CK



Complete Liner

10-15mmHg

Features

- Mild compression in foot and ankle with graduated compression in calf
- Comfortable knee-high liners
- *drirelease* moisture wicking and odor control
- Extra thick sole from toe to heel for added comfort
- Silver ionized treatment for bacteriostatic properties
- Available in Black and Beige

Transition Liner Sizing Chart (Included with Compreflex Complete)

SIZE	Ankle Circ. (A)	Calf Circ. (C)	Women's Shoe	Men's Shoe	Pair (Black)	Pair (Beige)
Small	19-26cm	29-39cm	<7½	<6½	1601-CK-F1	1611-CK-F1
Medium	21-30cm	34-44cm	7.5-9.5	6.5-8.5	1602-CK-F1	1612-CK-F1
Large	26-36cm	40-50cm	10-12	9-11	1603-CK-F1	1613-CK-F1
X-Large	31-41cm	46-56cm	12.5-15	11.5-14	1604-CK-F1	1614-CK-F1
XX-Large	36-46cm	54-64cm	12.5-15	11.5-14	1605-CK-F1	1615-CK-F1




Accutab Devices



Easy To Use

- Gently lay the Compreflex straps in place without stretching the fabric
- Attach one Accutab at the end of each strap
- Pull each strap to the prescribed dosage indicated by the lines on the Accutab: 20–30, 30–40, or 40–50mmHg



Compreflex Strap Extender

- Provides up to 10cm of additional circumferential size range to Compreflex calf wraps
- Detaches easily
- Potential solution for patients outside of standard size ranges of Compreflex calf wraps, or needing custom size

	1 Unit	10 Unit Pack
Black	1401-SE1	1401-SEB
Beige	1411-SE1	1411-SEB



Compreboot Lite

Variable compression



Features

- Ultra-thin profile that easily fits into most footwear
- Highly adjustable design (may be trimmed to fit)
- Ideal for mild edema/lymphedema
- Breathable stretch fabric conforms to ankle
- Provides mild, therapeutic compression
- Available in Black only



Compreboot Standard

Variable compression



Features

- 4-way stretch for firm compression around ankle
- Included with Compreflex Standard Calf & Foot
- Ideal for mild-moderate edema/lymphedema
- Features neoprene material with hook and loop tabs
- Open heel and toe
- Provides moderate, therapeutic compression
- Available in Beige and Black



Compreboot Plus

Variable compression




Features

- Chipped foam lining
- Ideal for chronic/acute edema/lymphedema
- Removable gripper patches provide traction
- Straps are removable for custom positioning or replacement
- Provides moderate, therapeutic compression
- Available in Black only



REQUIRED MEASUREMENTS



Measuring Instructions

I: Measure circumference of ankle bend and heel.

J: Measure circumference across metatarsal heads.

K: Measure length on the medial side of the foot from the heel to the 1st metatarsal head.

Compreboot Lite and Compreboot Standard Chart

	Small		Medium / Large		X Large / XX Large	
	Regular	Long	Regular	Long	Regular	Long
I	28-36cm	28-36cm	39cm max	39cm max	44cm max	44cm max
J	22-26cm	22-26cm	30cm max	30cm max	33cm max	33cm max
K	14-18cm	19-23cm	18-20cm	23-25cm	20-22cm	25-27cm

Compreboot Lite Product Codes

		Small/ Medium/ Large	X Large XX Large
Black	Regular	1803-BLR	1805-BLR
	Long	1803-BLL	1805-BLL

Compreboot Standard Product Codes

		Small	Medium Large	X Large XX Large
Black	Regular	1801-BTR	1803-BTR	1805-BTR
	Long	1801-BTL	1803-BTL	1805-BTL
Beige	Regular	1811-BTR	1813-BTR	1815-BTR
	Long	1811-BTL	1813-BTL	1815-BTL

Compreboot Plus Chart

	Small		Medium / Large		X Large / XX Large	
	Regular	Long	Regular	Long	Regular	Long
I	28-36cm	28-34cm	39cm max	39cm max	44cm max	44cm max
J	22-26cm	22-25cm	30cm max	28cm max	33cm max	33cm max
K	14-18cm	19-23cm	18-20cm	23-25cm	20-22cm	25-27cm

Compreboot Plus Product Codes

		Small	Medium/ Large	X Large XX Large
Black	Regular	1801-BPR	1803-BPR	1805-BPR
	Long	1801-BPL	1803-BPL	1805-BPL

GLOVES

Rubber Gloves* – Ridged
The practical and protective
way to apply and remove
your compression stockings.



Package: 12 per Pack

Size	Product Code
Small	591R400S
Medium	591R400M
Large	591R400L
X-Large	591R400X
ASST'D 4xS/M/L	591R400A

Textile Gloves* – These high-
quality gloves have a latex
surface to provide a good grip
when putting on and taking
off compression stockings.



Package: 12 per Pack

Size	Product Code
X-Small	593R400XS
Small	593R400S
Medium	593R400M
Large	593R400L
X-Large	593R400X
ASST'D 4xXS/S/M/L	593R400A

Glove Sizing Chart - Measure Circumference of the Hand at the Palm

16.5–18cm (6.5"–7")	19–20.5cm (7.5"–8")	21.5–23cm (8.5"–9")	24–25.5cm (9.5"–10")
Small	Medium	Large	X-Large



Latex-Free Gloves
Gripping donning gloves
recommended for use in a
latex-free environment.

Package: 12 per Pack

Size	Product Code
Small	592R400S
Medium	592R400M
Large	592R400L
X-Large	592R400X

Textile Glove Sizing Chart

Measure Circumference of the Hand at the Palm

Length	20.5–21cm	21.5–22cm	22.5–23cm	23.5–24cm	24.5–25cm
Width	7–7.5cm	7.5–8cm	8–8.5cm	9–9.5cm	10–10.5cm
X-Small	X-Small	Small	Medium	Large	X-Large

DONNING AID



Sim-Slide®

Size	Product Code
Small	587SLIDESM
Medium	587SLIDEMD
Large	587SLIDELG
X-Large	587SLIDEXL

Sizing Chart

For use with Open toe styles

NEW SIZING

Women's Shoe Size	≤4–6	6.5–8	8.5–9.5	≥10
Men's Shoe Size	3–4.5	5–6.5	7–8	≥8.5
Size	S	M	L	XL



Magnide On/Off®

Size	Product Code
Medium	587MAGMD
Large	587MAGLG
X-Large	587MAGXL

Sizing Chart - For use with Closed toe styles

Women's Shoe Size	≤8.5	9–12	≥12.5
Men's Shoe Size	≤7.5	8–11	≥11.5
Size	M	L	XL

* **Caution:** Rubber Gloves and Textile Gloves contain natural rubber latex that may cause allergic reaction



SIMON

A donning and doffing aid so versatile that it can easily be stowed to accompany you wherever you go.

Size	Calf Circumference	Product code
X-Small	25–29cm	589SX5
Small	30–38cm	589SSM
Medium	38–46cm	589SMD
Large**	46–54cm	589SLG



Washing Solution

Gentle washing solution specifically formulated to help maintain the quality of SIGVARIS compression socks and stockings. Also great for fine washables!

Package Options

Contains	Product Code
Six 16 Oz. Bottles (473ml)	586W900
One 16 oz. Bottle (473ml)	586W950
Twelve 2 oz. Bottles (60ml)	586W975



HOLD UP

Roll-on water-soluble liquid adhesive to keep stockings in place.

Package Options

Contains	Product code
Twelve 2 oz. Bottles (60ml)	5811400
One 2 oz. Bottle (60ml)	5811450



Easy Slide Arm

Application aid for armsleeves.

Size	Product Code
Small/Medium	587E01408
Large	587E01425



Cotton Underliner

10–15mmHg stocking. Can be used as an under layer to protect fragile skin, secure wound dressings, add 10–15mmHg of compression, or create an easy glide surface to aide in putting on and taking off compression stockings. For sizing, please refer to the Natural Rubber Calf sizing on page 43.

Package: 2 per pack

Size	Product Code
Small / Short	588USS
Small / Long	588USL
Medium / Short	588UMS
Medium / Long	588UML
Large-XL / Short	588ULS
Large-XL / Long	588ULL



SIGVARIS Doff N' Donner

A soft device that allows you to quickly guide a compression sock or stocking up and down your leg or arm, making it simple to put on and remove all styles of graduated compression. It is perfect for both self-donning and caregiver use.

Size	Product Code
One Size	583D100



Foot Slip

Comfortable applicator for use with open toe stockings.

Package: 6 per pack

Size	Product Code
One Size	580F200



SIGVARIS Doff N' Donner Cone

This loading device makes it easy to roll stockings onto the SIGVARIS Doff N' Donner. The bottom of the base firmly attaches to most clean, flat surfaces via suction.

Size	Product Code
One Size	583D200

Symptoms / Indications	Compression Level at the Ankle
<ul style="list-style-type: none"> Tired, achy legs Occupational or evening edema Leg discomfort from long hours of standing or sitting Prevent varicosities & swelling during pregnancy Preventative for legs predisposed to risk Reduce swelling during travel 	15–20mmHg
<ul style="list-style-type: none"> Heavy, fatigued, aching legs Mild edema in lower extremities Mild varicosities with minimal edema Mild varicosities during pregnancy Lymphedema upper extremities Post-procedure of small veins Prevention or management of DVT <p><i>Note: A prescription is recommended for stockings with the indications of 20mmHg or higher</i></p>	20–30mmHg
<ul style="list-style-type: none"> Moderate varicosities with mild to moderate edema (also during pregnancy) Post-procedure of larger veins to maintain treatment success CEAP C3 (edema), C4 (skin changes without ulceration), C5 (skin changes with healed ulceration), C6 (skin changes with active ulcerations) Management and treatment of DVT or Post-Thrombotic Syndrome Superficial phlebitis Lymphedema after decongestant therapy to maintain reduction Orthostatic / postural hypotension Post-fracture, post-traumatic edema 	30–40mmHg
<ul style="list-style-type: none"> Lymphedema after decongestant therapy to maintain reduction Severe Post-Thrombotic Syndrome Severe varicosities and/or edema 	40mmHg or higher
<ul style="list-style-type: none"> Severe arterial insufficiency, intermittent claudication, ischemia Uncontrolled congestive heart failure 	Contraindications
<ul style="list-style-type: none"> Skin sensitivities or allergies Advanced neuropathy with sensory loss in leg or foot Diabetes with arterial compromise Confinement to bed or non-ambulatory use unless otherwise prescribed by the physician Eczema, hypodermatitis, dermatosclerosis, dermatitis 	Cautions
<p>The therapeutic effect of graduated compression stockings is based on the compression level. With the proper CEAP classification, the correct compression level at the ankle can be chosen. For more severe symptoms, consult with a vascular specialist.</p>	
<p><i>No liability accepted for non-observance of contraindications and cautions.</i> Fronck HS, Bergan JJ, et al. The Fundamentals of Phlebology: Venous Disease for Clinicians. 2004. pg 151. Blaettner, et al. Leg Symptoms of Healthy People and Their Treatment with Compression Hosiery. Phlebology 2008. Compression therapy of the extremities. Partsch, Rabe, Stemmer. ISBN: 2.85480.770.7.</p>	

Term	Definition	Sigvaris Usage
CHRONIC VENOUS DISORDERS (CVD)	A disease state comprised of disorders ranging from tired, achy legs, varicose veins, skin issues to venous leg ulcers (open and healed).	SIGVARIS GROUP manufactures MEDICAL and CONSUMER HEALTH CARE products to help reduce the risk of, improve symptoms of, and provide treatment for CVD.
DOUBLE-COVERED YARN	Spandex fibres are covered twice with yarn to prevent bare spandex from touching the skin for improved wearing comfort.	Used in nearly all SIGVARIS GROUP graduated compression socks and hosiery to create products that are easy to put on and take off, while also providing greater comfort.
DEEP VEIN THROMBOSIS (DVT)	A blood clot in a deep vein, most commonly in the lower leg or thigh.	SIGVARIS GROUP cares about DVT prevention. Clinical evidence supports the use of graduated compression to help prevent DVTs and the more serious complication of PE (Pulmonary Embolism), which can be life-threatening.
GRADUATED COMPRESSION VS. ANTI-EMBOLISM	Graduated compression is for ambulatory patients, whereas anti-embolism stockings are for non-mobile patients.	SIGVARIS GROUP offers true graduated compression, meaning products are tightest at the ankle and decrease in pressure going up the leg to help improve circulation.
HYPOALLERGENIC	Designed to minimize the chance of an allergic reaction.	SIGVARIS GROUP offers latex-free products and cotton fabrics for those with sensitive skin and allergies. Also, many SIGVARIS GROUP accessories are designed for consumers with sensitive skin.
ISO CERTIFIED	SIGVARIS GROUP is ISO (International Organization for Standardization) certified, indicating we have the required quality systems in place.	SIGVARIS GROUP follows ISO Standard 13485:2003 for our products manufactured in North and South America.
MICROFIBER	An extremely fine synthetic fibre knitted into textiles with the texture and drape of natural-fibre fabric. This fabric has excellent breathability.	Used to create products with a soft and comfortable feel. Find this fibre in 840 Soft Opaque, 830, 143 & 183 Microfiber Patterns, 820 Microfiber, and 189 BUSINESS CASUAL.
mmHg (MILLIMETERS OF MERCURY)	Defined as the extra pressure generated by a column of mercury one millimeter high; the same measurement used for blood pressure.	SIGVARIS GROUP products are offered in several ranges of compression levels based on this measurement.
MERINO WOOL	Merino is the softest and finest wool available in the world and is inherently itch-free.	SIGVARIS GROUP uses this premium wool yarn to create soft products that are thermoregulating for use in all seasons. Used in 422 Thermoregulating Wool, 421 MERINO OUTDOOR PERFORMANCE SOCKS, 240 Merino Wool, and 152 & 192 ALL-SEASON MERINO WOOL.
SEA ISLAND COTTON	Sea Island Cotton is the rarest cotton in the world and represents the finest in cotton softness.	SIGVARIS GROUP uses the highest-quality yarns in the industry to create products that are best in class. Look for this fine quality cotton in our 220 Sea Island and 151 & 191 SEA ISLAND COTTON SOCKS.
SENSINNOV	A patented SIGVARIS silicone top band that gently grips the skin and keeps the band in place all day.	This top band is featured in 780 Sheer, 840 Soft Opaque, 710 Patterns, and 820 Microfiber to help keep thigh-high stockings in place all day long.
Supima COTTON SUPIMA WORLD'S FINEST COTTONS	Supima accounts for only about three percent of annual cotton production in the United States. Its fineness and longer staple length make it a premium cotton fibre.	SIGVARIS GROUP uses this premium cotton yarn in the 230 Cotton, 220 Sea Island and 151 & 191 SEA ISLAND COTTON SOCKS, 146 & 186 CASUAL COTTON, and 230, 142 & 182 Cushioned Cotton products, making them comfortable and ideal for those with sensitive skin.



**Sensitive
Skin**



Thermoregulating



Breathable



**Moisture
Wicking**



Soft



Durable



Cotton



Wool



Microfiber



**Natural
Rubber**



**Largest Size
Range**



**Wide Size
Range**



Sheer



**Medium
Sheer**



Opaque



Stylish



Stylish



Elegant



Colours



Patterns



**Sensinno
Grip-top
Band**



**Anatomically
Contoured**



**Cushioned
Sole**

Please follow these instructions carefully to determine the correct size for different product styles.

Measurements should always be taken in the following order, and as early in the morning as possible when swelling is minimal.

CALF STYLE

Open/Closed Toe

1	Measure the ankle circumference at the narrowest part of the ankle, above the ankle bone. Use the appropriate series sizing chart to determine if ankle is SMALL, MEDIUM, LARGE or X-LARGE. Note: This should be your first measurement.
2	Measure the calf circumference at the fullest part of the calf.
3	Determine the calf length by measuring from the floor up to the bend of the knee. This measurement is taken on the outside of the leg.

PANTYHOSE AND THIGH-HIGH STYLE

Open/Closed Toe

1	Measure the ankle circumference at the narrowest part of the ankle, above the ankle bone. Use the appropriate series sizing chart to determine if ankle is SMALL, MEDIUM or LARGE. Note: This should be your first measurement.
2	Measure the calf circumference at the fullest part of the calf.
3	Skip Step 3 (but refer to General Guidelines #6 below).
4	Measure the thigh at its widest circumference.
5	Determine the leg length by measuring from the floor up to the gluteal fold.
6	Determine the hip circumference by measuring the widest part of the hips.

GENERAL GUIDELINES

1	Measurements should always be taken as early in the morning as possible, when swelling is minimal.
2	Measure against bare skin, not over clothing (Steps 1–5).
3	The measuring system is based on having the foot flat against the floor and the calf or leg at a right angle to the foot (see images to right for examples).
4	Measurements should be taken by pulling the tape around the leg so that it is snug but not constricting.
5	Length measurements should be taken from the ground up.
6	Take all measurements, even if the patient only requires calf-length products, as the patient might decide at a later date to purchase additional styles.




1. Ankle circumference at narrowest part of ankle




2. Calf circumference at fullest part of the calf




3. Calf length from the floor to the right-angle bend of the knee



4. Widest circumference of the thigh



5. Leg length from the floor up to the gluteal fold



6. Circumference at the hip

How to Order Medical and Specialty Products

Series Number 23 = COTTON (See individual pages for series number)	Compression Level 1 = 15–20mmHg 2 = 20–30mmHg 3 = 30–40mmHg 4 = 40–50mmHg* 5 = 50–60mmHg*	Product Style C = Knee-high or Calf (See individual pages for available styles)	Size ML = Medium Long (See individual pages for sizing)	Gender M = Men W = Women O = Open toe	Colour 99 = Black (See individual pages for available colours)
23	2	C	ML	M	99
* 40–50mmHg and 50–60mmHg available in Natural Rubber only					
To order Plus Calf with Grip-top or Plus Size Pantyhose (860P for women only): add "/P" after the colour in the ordering code					
To order Calf Style with Grip-top: add "/S" after the colour in the ordering code (see product pages for availability)					
To order Thigh with Waist Attachment: add "/L" to indicate the Left leg or "/R" to indicate the Right leg after the colour in the ordering code (860W and 500W only)					
To order MOTION HIGH TECH and MOTION THERMOREGULATING WOOL: The letter for Gender is not used in the ordering code. Example: 412CSM59 = MOTION HIGH TECH, 20–30mmHg, Closed toe Calf, Size Small Medium, Colour Red					

How to Order CONSUMER HEALTH CARE Products

Series Number 120 = SHEER FASHION (See individual pages for series number)	Product Style N = Thigh-high (See individual pages for available styles)	Size A = Small (See individual pages for sizing)	Colour 33 = Natural (See individual pages for available colours)
120	N	A	33
To order ATHLETIC RECOVERY SOCKS, PERFORMANCE SLEEVES, and MERINO OUTDOOR SOCKS: Size is only one letter and the letter for Gender is not used in the ordering code. Example: 421CL31 = MERINO OUTDOOR SOCK, 15–20mmHg, Closed toe Calf, Size Large, Colour Olive			

To Order:
1.800.363-4999

Online Dealer Portal:
sigvaris-online.com

Email address:
canadaorders@sigvaris.com

WHY SIGVARIS GROUP?

We produce the best compression garments to enhance our customers' health and quality of life.

We are dedicated to helping legs look and feel their best. Using premium yarns, we've spun ourselves into being the worldwide leader in graduated compression socks and hosiery by embracing our core beliefs of QUALITY, COMFORT, FIT and INNOVATION. Staying true to these values, we are continuing our efforts by expanding our business to help customers with medically-complex needs.

More than 20 people handle our products from start to finish, ensuring every garment is of the highest quality.



We offer the most comfortable fabric options available, including the finest cottons, wools and microfibers. We double-cover the yarns for ease of donning and doffing.

We offer premium products designed for all shapes and sizes, from petite to plus, including compression for a broad spectrum of edemas. All of our products are designed using the most recent human data measurement studies. A proper fit is key to adherence.



We use the finest textiles and materials to create modern styles in our circular-knit products, while maintaining medical efficacy. We continue to invest in medical research to help customers.



The Worldwide Leader
in Compression Therapy

SIGVARIS, SIGVARIS GROUP, PFS², ACCUTAB, and Sensinnov are registered trademarks of SIGVARIS AG, CH-9014 St.Gallen/Switzerland, in many countries worldwide. © Copyright 2019 SIGVARIS, Inc.

HOLD UP is a trademark of Spec Enterprise LLC. Lycra is a registered trademark of INVISTA.

Supima is a registered trademark of the Supima Association of America Corporation. NILIT is a registered trademark of Nilit. All rights reserved.

drirelease is copyright ©2019 drirelease. drirelease is an Optimer Brands company. All Rights Reserved.

Sim-Slide and Magnide are registered trademarks of Arion Holding B.V.

Sigvaris Dealer Portal

SIGVARIS
GROUP

The SIGVARIS GROUP Dealer Portal is your online location for securely ordering SIGVARIS GROUP products and managing your account. Place and track your orders, review order history, contact Customer Care, and much more at sigvaris-online.com.

DEALER PORTAL BENEFITS


24/7 access – order any time!

Three ways to securely place an order online:

Catalog, Quick Entry or Order Upload

Access to order history and tracking

Email canadaorders@sigvaris.com
to get your account activated!



Sigvaris Education



SIGVARIS GROUP offers educational seminars across Canada in the Spring and Fall to help increase your leg health knowledge.

Become A SIGVARIS Certified Fitter!

Contact your Territory Manager

or call Customer Care at **1-800-363-4999**

or email canadacares@sigvaris.com

We are now on
social media!
Check us out
on Facebook
and Instagram!



sigvarisgroup.canada



[sigvarisgroup_canada](https://www.instagram.com/sigvarisgroup_canada)



Established in 1864, SIGVARIS GROUP is, and has always been, a family-owned company. Founded and headquartered in Switzerland, our passion for helping people feel their best is why we've been a global market leader for 155 years and counting.

Our compression garment portfolio includes MEDICAL, CONSUMER HEALTH CARE and products for MEDICALLY COMPLEX EXTREMITIES that are specifically designed to restore, assist and enhance people's lives. We are renowned for making exceptional products that combine advanced technology and quality craftsmanship with diverse styles that fit our modern customers' lives.

When heritage, competence and dedication come together, there's no limit to where we'll go.

Live your life. Feel your best.

For more information, please visit:
sigvaris.com

TO ORDER:

800.363.4999 (phone)
800.263.8736 (fax)
canadaorders@sigvaris.com

Online Dealer Portal:

sigvaris-online.com